


RESTORATION[®] MODULAR


stryker[®]

REVISION HIP SYSTEM

Surgical Protocol

**Restoration Modular
Calcar Body/Conical Distal Stem
Femoral Components
Using the Restoration Modular
Instrument System**





Restoration Modular
Revision Hip System

Restoration® Modular Revision Hip System

Surgical Protocol

**Restoration Modular
Calcar Body/Conical Distal Stem
Femoral Components
Using the Restoration Modular
Instrument System**

Indications

- Noninflammatory degenerative joint disease, including osteoarthritis and avascular necrosis;
- Rheumatoid arthritis;
- Correction of functional deformity;
- Revision procedures where other treatments or devices have failed; and,
- Nonunions, femoral neck fractures, and trochanteric fractures of the proximal femur with head involvement that are unmanageable using other techniques.
- The Restoration Modular Hip System is intended for primary and revision total hip arthroplasty as well as in the presence of severe proximal bone loss. These femoral stems are designed to be press fit into the proximal femur.

Contraindications

- Active infection or suspected latent infection in or about the hip joint;
- Bone stock that is inadequate for support or fixation of the prosthesis;
- Skeletal immaturity;
- Any mental or neuromuscular disorder that would create an unacceptable risk of prosthesis instability, prosthesis fixation failure, or complications in postoperative care.

See package insert for warnings, precautions, adverse effects, information for patients and other essential product information.

Before using Restoration Modular instrumentation, verify:

- Instruments have been properly disassembled prior to cleaning and sterilization;
- Instruments have been properly assembled post-sterilization;
- Instruments have maintained design integrity; and,
- Proper size configurations are available.

For instructions for Cleaning, Sterilization, Inspection and Maintenance of Orthopaedic Medical Devices, refer to LSTPI-B.

Acetabular Options

Stryker Orthopaedics offers a wide variety of acetabular components that are compatible with the Restoration Modular Femoral Components. The surgeon should refer to a specific acetabular component's surgical technique for a discussion of acetabular surgical procedures. The Restoration Modular Hip System is compatible only with Stryker Orthopaedics femoral bearing heads listed in the chart on page 16 and page 21.

Restoration® Modular

Surgical Protocol

System Overview

The Modular Calcar Body/Conical Distal Stem Femoral components are part of the Restoration Modular Revision Hip System. The system takes advantage of the long clinical experience with distally fixed implants, while making use of modern technology to enhance proximal load transfer to the femur. This is achieved by mating a selected proximal body with a selected distal stem to provide a femoral prosthesis that minimizes proximal-distal mismatching, often associated with monolithic implants.

Revision hip surgery is very complex in that the surgeon may face compromised soft tissues, retained cement, severe bone loss, and poor residual bone. A set of implant options is essential to best fit the implant to the present bone defect. The Restoration Modular Calcar Body/Conical Distal Stem Femoral components were designed specifically for use in revision cases in which the femoral bone stock is severely compromised in the proximal third or proximal half of the femur. They also may be used for less challenging reconstructive surgery ranging from difficult primary up to, and including, Type III revision cases.[†]


The titanium alloy (Ti-6Al-4V ELI) Calcar Bodies are circumferentially plasma sprayed with commercially pure titanium and then over-sprayed with PureFix HA. These surface enhancements have demonstrated biocompatibility through many years of use at Stryker Orthopaedics.* Proximally, the Calcar Body segment helps maintain rotational and axial stability when adjacent to viable bone. Seven Calcar Body diameters are available (range 19mm through 31mm in 2mm increments) with three vertical offsets: +10mm, +20mm, and +30mm. These vertical offsets may be used to adjust overall stem length intraoperatively. (See Sizing Charts on page 3.) Cobalt chrome bushings have been added to the vertical flanges of the Calcar Body to allow the utilization of 2.0mm Dall-Miles Beaded Cables for trochanteric reattachment and cerclage fixation.

The Calcar Bodies incorporate a V40 taper and are compatible with CoCr, Biolox delta Ceramic and Alumina Ceramic Femoral Heads.

The Conical Distal Stem is designed to provide immediate diaphyseal rotational and axial stability. These stems are available in three lengths – 155mm, 195mm, and 235mm. Each distal stem length comes in 15 fluted, conical diameters from 14mm to 28mm in 1mm increments. The 155mm and 195mm Conical Distal Stems are offered with a straight design option. The 195mm Conical Distal Stem is also offered with a bowed option. The 235mm Conical Distal Stem is only offered with a bowed option.

The total length of the Calcar Body/Conical Distal Stem construct will be dependent upon the body and stem chosen. Stem lengths are measured from the chosen Calcar Body (+10mm, +20mm, +30mm) with a +0mm (STD) Femoral Head from the head center to the distal tip of each of the three lengths of Conical Distal Stems. Review Sizing Charts on page 3 for stem lengths.

Note: The Calcar Body/Conical Distal Stem lengths are measured using the +10mm Calcar Body with a +0mm (STD) Femoral Head from the head center to the distal tip of the Conical Distal Stem.



[†] D'Antonio, J., et al. Classification of Femoral Abnormalities in Total Hip Arthroplasty. *Clin Ortho and Rel Research*. 1993; Number 296: pp. 133 – 139. Longjohn, D. & Dorr, L. Bone Stock Loss and Allografting: Femur. *Revision Total Hip Arthroplasty*. 1999. pp. 100 – 111.

* Capello WN, D'Antonio JA, Feinberg JR, Hydroxyapatite in Total Hip Arthroplasty: Clinical Results and Critical Issues. *CORR*. 355: 200-211. 1998.

Restoration® Modular

Surgical Protocol

Stem Length Options

Calcar Body Sizes and Head Offsets with V40 Femoral Heads
available in 22mm, 26mm, 28mm, 32mm, 36mm, 40mm & 44mm

Calcar Body Sizes	-4mm*	+0mm (STD)	+4mm	+8mm	+12mm
19mm	31mm	34mm	37mm	40mm	43mm
21mm	33mm	36mm	39mm	42mm	45mm
23mm	37mm	40mm	43mm	46mm	49mm
25mm	41mm	44mm	47mm	50mm	53mm
27mm	41mm	44mm	47mm	50mm	53mm
29mm	41mm	44mm	47mm	50mm	53mm
31mm	41mm	44mm	47mm	50mm	53mm

*Not available in 22mm or 26mm diameter head (see Head Compatibility chart on pages 16 or 21).

Calcar Body/Conical Distal Stem


Calcar Body Sizes	Neck Angle	Distal Stem Lengths (mm)	Distal Stem Diameters
19mm 21mm 23mm 25mm 27mm 29mm 31mm	132°	155, 195, 235	14mm – 28mm in 1mm Increments

Calcar Body/Conical Distal Stem Combined Overall Lengths**

Conical Distal Stem Sizes	Calcar Body Heights		
	80mm +10mm	90mm +20mm	100mm +30mm
155mm Length 14mm-28mm dia. (Straight)	235mm	245mm	255mm
195mm Length 14mm-28mm dia. (Straight & Bowed)	275mm	285mm	295mm
235mm Length 14mm-28mm dia. (Bowed)	315mm	325mm	335mm

**Femoral head neck length options will increase overall stem lengths – range -4mm, +0mm (STD), +4mm, +8mm, and +12mm.
Head center (+0mm STD) to distal stem tip.

Cobalt Chrome Bushings are for use with 2.0mm Dall-Miles Beaded Cables.



Restoration® Modular

Surgical Protocol

Bone Defect Classifications

Type 1 - Minor Bone Loss

- The metaphysis is expanded, but intact.
- The calcar is partially absent.
- There is minimal bone loss anteriorly and posteriorly.
- The diaphysis is intact.

Type 2 - Significant Bone Loss

- The metaphysis is compromised.
- There is no calcar.
- There is minimal bone loss anteriorly and posteriorly.
- The available proximal bone may be thin, sclerotic, and incapable of support.
- The diaphysis is intact.

Type 2A - The calcar is non-supportive, but the diaphysis is still intact.

Type 2B - The calcar is non-supportive, the anterolateral metaphysis is deficient, but the diaphysis is still intact.





Type 2C - The calcar is non-supportive and the posteromedial part of the metaphysis is deficient, but the diaphysis is still intact.

Type 3 - Massive Bone Loss

- Complete circumferential bone loss in the metaphysis, extending to the diaphysis.
- The metaphysis and part of the diaphysis are deficient.
- The anterolateral bone and supporting subtrochanteric metaphyseal bone are absent.
- The metaphysis is not stable and will not offer rotational stability.
- There is massive bone loss anteriorly and posteriorly.
- The stability of the implant is dependent on distal diaphyseal fixation.

Type 4 - Massive Bone Loss


- Extensive circumferential segmental bone loss proximally.
- Extensive cavitory loss involving the entire diaphysis.
- Extensive ectasia of the diaphysis.
- Proximal femoral allograft required with reduction osteotomy of the diaphysis.
- Cortical diaphyseal bone is often thin and needs to be supplemented with cortical strut grafts.
- Segmental defects can be repaired with cortical strut graft and cerclage wiring, and cavitory defects can be filled with impacted particulate graft.



Source: D'Antonio, J., et al. Classification of Femoral Abnormalities in Total Hip Arthroplasty. *Clin Ortho and Rel Research*. 1993; Number 296: pp. 133 – 139.
 Longjohn, D. & Dorr, L. Bone Stock Loss and Allografting: Femur. *Revision Total Hip Arthroplasty*. 1999. pp. 100 – 111.

Restoration® Modular

Surgical Protocol



Preoperative Evaluation and Planning

The Restoration Modular Calcar Body/Conical Distal Stem Femoral Hip System offers a complete set of femoral X-ray templates for the surgeon to help assess the implant requirements. All seven Calcar Body Templates (with three vertical offsets each) can be combined with each of the Conical Distal Stem Templates. All templates are at 120% magnification. The use of mag markers will facilitate accurate magnification measurements. If mag markers are not used, measure the existing implants on the X-ray to ensure that magnification is approximately 120%.



Preoperative planning is strongly recommended for leg length planning, measuring the length of the existing prosthesis being revised, predicting the potential use and type of trochanteric osteotomy, the Calcar Body size and vertical offset, and the Conical Distal Stem diameter and length of the prosthesis to be implanted.

Anterior-Posterior (A/P) and Medial-Lateral (M/L) radiographs are necessary for X-ray templating. In cases of severe femoral compromise, a full A/P pelvic X-ray of the operative side as well as the contralateral side is helpful to assess the biomechanical requirements of the reconstruction. The lateral X-ray is informative in that it will show the anterior bow of the femur, which is useful when templating with the 155mm straight, 195mm straight and bowed, or 235mm bowed long stems.

First, position an acetabular template over the A/P radiograph, aligning the acetabular shell surface with the subchondral bone. Mark the center of rotation of the acetabulum indicated on the template.

Place the appropriate two-piece femoral template on the radiograph. Ensure that the distal length of the prosthesis will be sufficiently anchored in good cortical bone – this is generally two-to-three canal diameters below the tip of the existing implant or defect. The necessary proximal body height is selected to anatomically correct the leg length.

IMPORTANT: Adjunctive proximal fixation/support is required for stem diameters of 11mm, 12mm, and 13mm, and is recommended for stem diameters of 14mm and larger.



Restoration® Modular

Surgical Protocol

Patient Selection

Proper implant selection is critical to the stability and longevity of the femoral stem implant in hip arthroplasty. Proper implant selection must consider design, fixation, and environmental variables including: patient weight, age, bone quality and size, activity level and preoperative level of health, as well as the surgeon's experience and familiarity with the implant device. Longevity and stability of the implant may be affected by these factors. Surgeons should advise patients of these factors.

The smaller sized femoral stem implants are intended for use in patients with smaller intramedullary femoral canals. Their geometry has been reduced to accommodate the anatomy of the smaller intramedullary femoral canal, which thereby decreases their fatigue-strength and load-bearing characteristics. Therefore, patients with high physical activity levels, poor bone quality, or who are overweight may be poor candidates for the smaller femoral implant stem.

Patients with high-activity level and/or higher weight patients are at greater risk for implant complications or failures. For patients with poor proximal bone quality, the use of supplemental adjunctive proximal fixation/support is advised for implant stability.


The surgeon must evaluate each situation carefully based upon the patient's clinical presentation before making any decisions regarding the selection of the implant.



A full range of implant sizes provides choice in selecting an implant to meet the specific demands of each patient.

Restoration® Modular

Surgical Protocol



Determine the Approximate Implant Size

Note that the tip of the greater trochanter is approximately at the same level as the center of rotation of the femoral head.

If no change in leg length is necessary, then the Calcar Body and Femoral Head center that is closest to the center of rotation marks the appropriate neck length and femoral head offset required. If leg lengthening is required, choose the Calcar Body height, neck length and offset that places the center of the femoral head on the overlay above the center of rotation. If it is necessary to shorten the length of the femoral neck, then select the Femoral Head center below the center of rotation.

Once the proximal geometry has been determined, select the appropriate Conical Distal Stem diameter of the implant by establishing the region of the femoral cortices that appears to be perfectly defined or free from defects that will allow the implant to achieve 10cm - 12cm of suitable distal fixation. Determine also the length required to place the distal stem tip two-to-three canal diameters below the lowest distal defect.

IMPORTANT: Do not plan to use the +30mm Calcar Body or the +12mm Femoral Head preoperatively. Use the next larger diameter Conical Distal Stem in the same implant length so that additional vertical offset, neck length, and femoral head offset options are available for adjusting leg length intraoperatively.

Patient Positioning and Surgical Approach

Revision total hip surgery presents challenges not seen in primary surgery. Therefore, each surgeon should position the patient and use the surgical approach for revision total hip arthroplasty with which he is most familiar. Patient positioning, prepping and draping, the skin incision, soft tissue dissection, and hip dislocation are performed according to the surgeon's preferred technique, making certain to adequately expose the acetabulum and femur as required by each revision situation.

There are also many femoral and trochanteric osteotomy techniques available to surgeons that assist in implant removal, overall reconstruction, and finally, postoperative management. The surgeon should use osteotomies that he is most familiar with and that best fit the challenge faced by each particular revision situation.

Note: To reduce the potential for femoral fracture, it is recommended that areas of defects in the femur are prophylactically cabled prior to reaming and stem insertion. Dall-Miles Beaded Cables work well to assist the surgeon in this step.



Restoration® Modular

Surgical Protocol

Cement Removal

Implant removal and subsequent cement removal can be a challenging proposition. Surgeons should utilize methods they are most familiar with or are most appropriate for the many revision situations that may arise. The Gray Revision Instruments are helpful in removing the existing acetabular and femoral prostheses as well as bone cement if present.

After removal of the femoral component, the acetabular component is removed and the acetabulum is prepared. Cement and fibrous tissue still present in the femoral canal may be left to help minimize blood loss during acetabular preparation. After the acetabulum has been prepared, any remaining cement, fibrous tissue, or debris in the femoral canal may be removed and reaming begun.



Gray Revision Instruments

Calcar Neck Resection Guide

A Calcar Neck Resection Guide is available for those instances where a surgeon prefers to make the calcar resection prior to distal femoral preparation. This method can be utilized in a primary surgery, or to excise additional bone in a revision scenario (Figure 1). Alternatively, the calcar resection can be accomplished after final distal reaming, using the Post Resection Cutting Guide (see related section on page 15).

The resection level should be identical to the level chosen during preoperative templating. Key features of the Neck Resection Guide (Figure 2):

1. The horizontal groove in the upper portion of the slotted area helps to reference the proximal tip of the greater trochanter. This is a good landmark that generally coincides with the center of rotation for the femoral head. Align the horizontal groove with the tip of the trochanter. The notches on the medial finger of the guide correspond with the head centers of the noted diameters.
2. The horizontal resection level slots provide a plane for marking the level of the cut, or can be used as a cutting surface for the saw blade. Determine the appropriate neck resection level (+10mm, +20mm, or +30mm).
3. The long tail of the guide is used for alignment with the femoral shaft axis. It is designed to be inserted under the soft tissues of the posterior aspect of the femur.
4. Mark the desired resection level by passing a sagittal saw blade through the resection guide. Remove the guide and complete the cuts ensuring that the proper cut angle is maintained.

Note: The total construct length including proximal Calcar Body and distal stem must be determined *PRIOR* to femoral preparation.

Figure 1






Figure 2



Restoration® Modular

Surgical Protocol



Box Chisel and Starter Awl

The Box Chisel may be used to open the proximal femur prior to use of the Starter Awl or in conjunction with the Starter Awl.

Box Chisel Use Prior to the Starter Awl

After the osteotomy has been performed, the Box Chisel is introduced into the anatomic axis of the femoral shaft (Figure 3). This will remove a wedge of bone at the medial base of the greater trochanter, helping to achieve neutral/lateral alignment of the Starter Awl.

Use of the Starter Awl and Depth Markings

The Starter Awl can be used by hand or on power. It is designed to open the femoral canal to a diameter of 9.5mm. Assemble the T-Handle or Power Reamer to the proximal end of the awl and target the piriformis fossa to open the canal. The awl is very sharp; therefore, care must be taken to centralize the awl within the femoral canal before reaming is started, avoiding extra osseous penetration with the tip (Figure 4).

There are two bold depth marking grooves on the Starter Awl (200mm and 240mm), and markings for the +10mm, +20mm, and +30mm resection levels. Measurement for depth insertion of the Starter Awl when used with all Calcar Body/Conical Distal Stems is at the tip of the greater trochanter.

Note: Since there is no +0mm (STD) Calcar Body, all reaming must reach a minimum depth of 210mm (200mm + 10mm on the reamer) or 250mm (240mm + 10mm on the reamer).

Figure 3





Figure 4



Restoration® Modular

Surgical Protocol

Box Chisel and Starter Awl (continued)

Box Chisel Use With the Starter Awl

After the awl has been used to open the femoral canal, the T-Handle or Power Reamer is removed with the awl engaged in the isthmus of the femoral canal. The shaft of the awl may now be used as an axial guide coinciding with the long axis of the femur. The Box Chisel is cannulated so that it slides over the shaft of the awl, removing a wedge of bone at the medial base of the greater trochanter (Figure 5).

Reaming with the Conical Distal Reamers progresses sequentially after use of the Starter Awl.

Note: To reduce the potential for femoral fracture, it is recommended that areas of defects in the femur are prophylactically cabled prior to reaming and stem insertion. Dall-Miles Beaded Cables work well to assist the surgeon in this step.

Clear Out Reamer Use

The Clear Out Reamer is used to open up the proximal portion of the canal when preparing the 14mm Conical Distal Stems (both straight and bowed). The Clear Out Reamer is used after the Starter Awl and before the Conical Distal Reamers (Figure 6). The function of this reamer is to prepare the canal to accept the initial 19mm Proximal Cone Reamer.

The reamer is inserted into the canal until the correct depth marking on the shaft aligns with the tip of the greater trochanter. When preparing for the Calcar Body, the line corresponding to the preoperatively templated Calcar Body (+10mm, +20mm, or +30mm) should be aligned with the tip of the greater trochanter.

IMPORTANT: Adjunctive proximal fixation/support is required for stem diameters of 11mm, 12mm, and 13mm, and is recommended for stem diameters of 14mm and larger.

Figure 5






Figure 6





Distal Reaming

Use of the Conical Distal Reamer – 155mm, 195mm, 235mm Stems

Conical distal reaming for the 155mm, 195mm, or 235mm Conical Distal Stems can be accomplished by use of a T-Handle (Figure 7) or on power (Figure 8). Select the diameter of a Conical Distal Reamer starting with a size one or two millimeters smaller than the templated size. The reamer diameters are available in 1mm increments from 13mm - 28mm. There are three depth marking grooves on the shaft of the Conical Distal Reamers (225mm, 265mm, 305mm) which correspond to the distance from the tip of the greater trochanter to the tip of the 155mm Conical Distal Stem, 195mm Conical Distal Stem, or 235mm Conical Distal Stem, respectively (Figures 9 and 10).

Note that the tip of the greater trochanter is approximately at the same level as the center of rotation of the femoral head. Since the depth markings correspond to the distance from the center of a +0mm (STD) Femoral Head on a +0mm (STD) body, you must measure and mark the desired Calcar Body height +10mm, +20mm, or +30mm on the reamer as there is no +0mm (STD) Calcar Body (Figures 7 and 8).

If the greater trochanter is off or not present, the measurements made during preoperative templating are necessary to determine the approximate location of the greater trochanter or head center. Alternately, measurements may be taken from an X-ray of the contralateral side.

Note: For the 155mm Conical Distal Stems, ream to 225mm +10mm, +20mm, or +30mm, as it corresponds to the desired Calcar Body height.

Note: For the 195mm Conical Distal Stems, ream to 265mm +10mm, +20mm, or +30mm, as it corresponds to the desired Calcar Body height.

Note: For the 235mm Conical Distal Stems, ream to 305mm +10mm, +20mm, or +30mm, as it corresponds to the desired Calcar Body height.

Note: If you have reamed too deep, back out the reamer with the T-handle and resect to the desired predetermined level.

Figure 7




Figure 8




Figure 9




Figure 10



Restoration® Modular

Surgical Protocol

Distal Reaming (continued)/Implant Insertion

Use of the Conical Distal Reamer – 155mm, 195mm, 235mm Stems (continued)

Progressively ream until resistance accompanied by cortical chatter is encountered. The reamers must be advanced into the femoral canal until the appropriate depth aligns with the tip of the greater trochanter, or approximate center of rotation (Figure 11). It is important not to over-insert the Conical Distal Reamers as these are matched to a specific sized distal stem.

If good cortical contact is not achieved, increase the reamer diameter in 1mm increments and insert only as deep as required based on desired Calcar Body height.

SUGGESTION: As in the Preoperative Planning Section, it is required that the femoral canal be reamed to the predetermined total construct length.

Note: Intraoperative X-rays are valuable to gauge the position of the Conical Distal Reamers relative to the A/P and M/L femoral cortices and to the anterior bow of the femur.

Implant Insertion – Distal Stem

Thread the appropriate Conical Distal Stem onto the Distal Stem Inserter. The distal end of the inserter has a hex geometry with a spring-loaded threaded end that mates with a corresponding geometry on the stem. Make sure that the distal tip of the Distal Stem Inserter is correctly aligned with the hex orientation feature of the insertion hole of the implant (Figure 12). Fully and securely attach the instrument to the distal stem by turning the locking knob clockwise.

Figure 11






Figure 12



Restoration® Modular

Surgical Protocol



Implant Insertion – Distal Stem

Insert the Distal Stem

There are two options for inserting distal stems, the Version Control Stem Inserter (Figure 13A) and the Distal Stem Inserter (Figure 13B). Both inserters have depth groove markings that correspond to the center of a +0mm (STD) Femoral Head implant on each of the seven Calcar Bodies (+10mm, +20mm, and +30mm) (Figure 13C). The top three markings on the inserter correspond with the +10mm, +20mm, and +30mm resection levels.

Important: Note the location of the +10mm mark on the stem inserter.

Note: Preoperative planning should have ensured that the tip of the distal stem will pass any distal defects by two to three canal diameters and will have 10cm - 12cm of satisfactory mechanical stability. Make sure that sufficient distal fixation is attained with all Conical Distal Stems, especially those that are significantly larger than the templated stem size.

Impact the Conical Distal Stem into the femoral canal until the stem will not advance any further, achieving rotational stability and preventing subsidence. View the depth grooves on the stem inserter. These will align with the tip of the greater trochanter and must correspond with the predetermined body length that will be used (+10mm, +20mm, or +30mm).

The Conical Distal Stem may advance further into the canal than originally templated and reamed. If this occurs, one of the following techniques may be suitable:

- Less than 10mm of advancement can be made up with Femoral Head Trials and implants.
- Impact the stem to the next calcar resection level and reattach Calcar Guidepost and resect using the Post Resection Guide.
- Greater than 30mm of advancement, it is recommended to remove the distal stem and ream to the next larger diameter that will achieve rotational (and axial) stability.

Remove the stem inserter by turning the locking knob counter-clockwise on top of the inserter.

Note: Depending on the bow of the femur, the trunnion of a *Straight Conical Distal Stem* may sit against the anterior femur upon insertion.

The *Bowed Conical Distal Stem* is designed to move the trunnion off the anterior cortex in the same type of femur (Figure 14).

Figure 13A




Figure 13B




Figure 13C





Figure 14



Restoration® Modular

Surgical Protocol

Calcar Body Preparation

The Calcar Bodies are prepared by Proximal Cone Reamers which are available in 7 diameters: 19mm, 21mm, 23mm, 25mm, 27mm, 29mm, and 31mm. These diameters are measured at the most proximal level of the coating on the medial side of the Calcar Body implant.

Insertion of Proximal Calcar Reamer Guidepost

Remove the threaded Proximal Calcar Reamer Guidepost from the tray and thread it into the top of the implanted Conical Distal Stem until fully seated. Use the 5mm Hex Driver and the small or large T-Handle to ensure full seating of the Guidepost on the distal stem; excessive torque is not required when tightening (Figure 15).

Insertion of Proximal Cone Reaming

Starting with the 19mm Proximal Cone Reamer, commence proximal preparation for the Calcar Body.

Attach the Proximal Cone Reamer to a power source and advance it over the Proximal Calcar Reamer Guidepost until it bottoms out on the post and it is impossible to advance the reamer further – visualize this by looking at the alignment groove (on the Guidepost and reamer) or view the top of the proximal slot, which when fully seated, will show no gap between the reamer and the Guidepost (Figure 16).

Proximal Cone Reaming progresses in 2mm increments until satisfactory contact within the trochanteric region is felt. Make note of the +10mm, +20mm, and +30mm reamer grooves, using the tip of the greater trochanter or other landmark, since this will dictate the Calcar Body height to be used (Figure 17).

Generally the depth groove of the Cone Reamer corresponds to the measurement taken during preoperative templating, however, it may be one of the other levels.

Note: If a calcar resection was made using the Calcar Resection Guide, the depth groove of the Cone Reamer should correspond to the calcar resection cut made during femoral preparation; i.e., a +10mm calcar resection must correspond to a +10mm Proximal Cone Reamer depth marking (refer to Neck Resection, page 8).

WARNING: Failure to fully seat the Proximal Calcar Reamer Guidepost, or failure to fully bottom out the Proximal Cone Reamer on the Guidepost may prevent proper preparation of the bone for the Calcar Body.

Figure 15




Figure 16





Figure 17



Restoration® Modular

Surgical Protocol



Calcar Post Resection Cutting Guide

Insert the Conical Distal Stem to its final seating level (refer to Insert the Distal Stem, page 13). Attach the Calcar Resection Guidepost.

IMPORTANT: Note the final seating level of the implanted distal stem as this will determine the calcar resection level. Ideally, this should correspond to the preoperatively templated construct.

Two Post Resection Guides allow for left/right hip anatomy and surgical approach. Choose the Post Resection Guide so that the medial side of the A/P adjuster arm is positioned medially (Figure 18).

Slide the appropriate Post Resection Guide over the Calcar Resection Guidepost until an audible and tactile click at the most distal groove is heard and felt (Figure 19). Slide the cutting block over the A/P adjuster bar into one of the grooves that is appropriate for the resection (Figure 20).

Adjust the version of the Calcar Resection Assembly on the Guidepost. Tighten the ring adjuster firmly with the thumbscrew (Figure 20).

Use the sagittal saw to cut through the vertical and horizontal slots (Figure 20). Manually control the version of the Post Resection Guide to achieve the desired resection. Ensure the vertical cut clears the top of the greater trochanter.

Note: The saw blade may come in contact with the Guidepost or Post Resection Assembly.

Remove the Post Resection Assembly from the femur, and complete the vertical and horizontal cuts through the femur.

Note: A complete cut through the femur cannot be accomplished until the Guidepost is removed.

Figure 18




Figure 19





Figure 20



Restoration® Modular

Surgical Protocol

Calcar Body Trial

Assemble the Appropriate Calcar Body Trial to Conical Distal Stem

Select the Calcar Body Trial corresponding to the previously made resection cut and final Proximal Cone Reamer diameter. Assemble the 8mm Hex Locking Bolt Driver Shaft to the T-Handle (Figure 21). Position the appropriate Calcar Body Trial with the integral locking bolt onto the distal stem. Determine the appropriate version for the trial and then tighten the locking bolt with the locking bolt screwdriver assembly or the Distal Stem Insertter (Figure 22). Excessive torque is not required when tightening.

Attach Head Trial

Select the head diameter according to surgeon preference. The Femoral Head Trials have a circumferential groove, which identifies the level of the center of rotation. Select the appropriate Femoral Head Trial based on preoperative templating from the chart below. Attach the Femoral Head Trial to the Calcar Body Trial (Figure 23).

Head	Head Size	Head Offsets
CoCr V40	22	+0, +3, +8
	26	-3, +0, +4, +8, +12
	28	-4, +0, +4, +8, +12
	32	-4, +0, +4, +8, +12
	36	-5, +0, +5, +10
	40	-4, +0, +4, +8, +12
	44	-4, +0, +4, +8, +12
Alumina V40	28	-2.7, +0, +4
	32	-4, +0, +4
	36	-5, +0, +5
Alumina C-Taper (when used with C-Taper Adaptor Sleeve-catalog #17-0000E)	28	-2.5, +0, +5
	32	-2.5, +0, +5
	36	-5, +0, +5
delta BIOLOX V40	28	-4, -2.7, +0, +4
	32	-4, +0, +4
	36	-5, -2.5, +0, +2.5, +5, +7.5
delta Biolox C-Taper (when used with C-Taper Sleeve - catalog #17-0000E)	28	-2.5, +0, +2.5, +5
	32	-2.5, +0, +2.5, +5
	36	-5, -2.5, +0, +2.5, +5, +7.5
delta Biolox Universal Taper (when used with Universal Taper Sleeve - catalog #6519-T-XXX)	28	-2.5, +0, +4
	32	-2.5, +0, +4
	36	-2.5, +0, +4
	40	-2.5, +0, +4
	44	-2.5, +0, +4

Figure 21

Figure 22

Figure 23

Figure 24

Perform a trial reduction and assess the hip for stability, leg length, and overall range of motion. Adjust the Calcar Body Trial as necessary to achieve maximum joint stability. Mark the desired anteversion on the bone with methylene blue or with a Bovie, in line with the neck. Carefully remove the Head Trial and Calcar Body Trial.


If additional leg length is required, a longer proximal body may be used so long as the current trial is not a +30mm length body.

Additionally, if the trial reduction indicates that a Femoral Head with a "skirt" is required, it may be possible to increase the body height by 10mm and use a shorter Femoral Head to produce an equivalent neck length without using a "skirted" Femoral Head. This may be beneficial in increasing range of motion.

Note: If the Calcar Body Trial becomes fixed in the canal, it may be removed with a bone hook or Trial Body Removal Device (Figure 24).

Restoration® Modular

Surgical Protocol



Calcar Body Insertion/Taper Lock Gauge

Calcar Body Impaction

Based on the trial reduction, select the appropriate height Calcar Body implant. Lavage the area surrounding the proximal taper of the distal stem. Wipe the Conical Distal Stem trunnion clean, and align the neck and trunnion of the Calcar Body implant with the methylene blue marking or Bovie mark, indicating the desired anteversion on the distal stem trunnion. Attach the Proximal Body Impactor to the Calcar Body (Figure 25) and impact the Calcar Body implant onto the trunnion of the Conical Distal Stem maintaining proper anteversion (Figure 26). The impaction of the Calcar Body onto the trunnion of the Conical Distal Stem cold-welds the tapers, locking the components together.

Taper Lock Gauge

After the Calcar Body has been impacted onto the distal stem, the Taper Lock Gauge can be used to assess proper engagement of the body with the stem. Insert the Taper Lock Gauge through the proximal body until it is seated on the distal stem (Figure 27). Slide the handle down until it is fully seated in the proximal body (Figure 28). The slotted indicator on the top of the handle will align within the groove corresponding to the Calcar Body height implanted (+10mm, +20mm, +30mm) (Figure 29).

Note: If the indicator is outside the corresponding groove, it may be necessary to further impact the body, or re-ream the proximal femur to clear out any bone stock that may interfere with the body properly seating on the stem.

Figure 25




Figure 26




Figure 27




Figure 29





Figure 28



Restoration® Modular

Surgical Protocol

Locking Bolt Assembly and Tightening

Locking Bolt Assembly and Tightening

Place the Locking Bolt into the Calcar Body and tighten the Locking Bolt with the 5mm Hex Locking Bolt Driver assembly (Figure 30). Assemble the Torque Wrench and Torque Wrench Adapter, and apply a minimum of 150in-lb and a maximum load of 180in-lb torque to ensure that the Locking Bolt is sufficiently tightened (Figure 31). The Calcar Body Steady Handle must be used to hold the anteversion of the Calcar Body in place while applying torque. The steady handle counter balances the torque applied to the bolt to ensure that only the implant and not the femur is torqued.

IMPORTANT: For the Calcar Bodies, when the body and stem tapers are fully engaged, the entire head of the locking bolt will be seen 1mm-2mm above the shoulder of the Calcar Body (Figure 32).

Note: The Conical Distal Stems have Spirallock threads that will not loosen if the Locking Bolt is sufficiently tightened. The Spirallock thread form reduces vibration loosening, provides a more uniform load distribution, reduces stress concentration, reduces fatigue failure, and eliminates the need for additional locking devices such as end caps.

Bone Grafting

Femoral deficiencies should be planned for and appropriately addressed as discussed in the preoperative planning part of this protocol. If the femoral cortex above the diaphyseal stem fixation point is deficient, the surgeon should be prepared to apply cortical strut grafts to repair and strengthen the femur.

Figure 30




Figure 31






Figure 32



Restoration® Modular

Surgical Protocol



Preparing for Cable Use

Preparing for Cable Use

Each Calcar Body has holes fitted with Cobalt-Chrome bushings that allow the use of Dall-Miles 2.0mm Beaded Cables to pass through the vertical flanges of the body. The +10mm Calcar Body has four available holes, and the +20mm and +30mm Calcar Bodies have six. The Dall-Miles Recon and Trauma Cable System provides the surgeon with a means of achieving trochanteric reattachment and a variety of methods for cerclage fixation.

The holes are angled outward to allow for cabling to pass around the greater trochanter. Additionally, drill guides may be used to prepare the bone to allow cables to pass through the bone.

Drilling may be performed either medially using the Calcar Drill Guide, or laterally using the Lateral Drill Guide.

IMPORTANT: Do not use Non-Beaded Cables with the Calcar Body. A Beaded Cable is required for each hole used. Cables should NEVER be wrapped around the implant body or neck.

Method 1 - Medial Preparation

Preparation from the medial to the lateral side is performed using the Calcar Drill Guide and a 2.2mm drill. Place the tip of the 2.2mm drill into the Calcar Drill Guide (Figure 33). Be sure to keep the drill away from the neck of the implant. Insert the drill guide until the tip seats in the hole (Figure 34). Drill through the holes and completely through the bone on the lateral side of the implant, ensuring that the proper drilling angle is achieved (Figure 35). Pass the non-beaded side of the cable through the hole from the medial side of the implant to the lateral side of the bone.

Figure 33




Figure 34





Figure 35



Restoration® Modular

Surgical Protocol

Preparing for Cable Use (continued)

Method 2 - Lateral Preparation

Preparation from the lateral side to the medial side is performed using the Lateral Drill Guide. Turn the knob counter-clockwise to open the Lateral Guide Arm and allow for placement over the greater trochanter. Insert the Centering Post on the fixed portion of the guide into the holes on the medial side of the implant. Place the Lateral Drill Guide into the proper orientation and turn the knob clockwise to secure the guide on the bone. Use the 3.2mm drill and drill through the hole in the Lateral Guide Arm up to the implant, ensuring the proper drilling angle is achieved (Figure 36). Pass the non-beaded side of the cable through the hole from the medial side of the implant to the lateral side of the bone (Figures 36A and 36B).

See the Dall-Miles surgical technique for details on the use of this system.

Note: It is important to secure the bone on the lateral side of the implant when drilling to facilitate the passage of the cables.

IMPORTANT: Do not use Non-Beaded Cables with the Calcar Body. A Beaded Cable is required for each hole used. Cables should **NEVER** be wrapped around the implant body or neck.

Figure 36




Figure 36A






Figure 36B





Final Trial Reduction

Attach Head Trial

Select the head diameter according to surgeon preference. The Head Trials have a circumferential groove, which identifies the level of the center of rotation (Figure 37). Select the desired Femoral Head Trial based on trial reduction from the chart below. Attach the Head Trial to the Calcar Body. The head center of the Head Trial, when attached to the implant construct, should correspond with the tip of the greater trochanter.

At this point, a final trial reduction can be performed using the attached Femoral Head Trial.

Head	Head Size	Head Offsets
CoCr V40	22	+0, +3, +8
	26	-3, +0, +4, +8, +12
	28	-4, +0, +4, +8, +12
	32	-4, +0, +4, +8, +12
	36	-5, +0, +5, +10
	40	-4, +0, +4, +8, +12
	44	-4, +0, +4, +8, +12
Alumina V40	28	-2.7, +0, +4
	32	-4, +0, +4
	36	-5, +0, +5
Alumina C-Taper (when used with C-Taper Adaptor Sleeve-catalog #17-0000E)	28	-2.5, +0, +5
	32	-2.5, +0, +5
	36	-5, +0, +5
delta BIOLOX V40	28	-4, -2.7, +0, +4
	32	-4, +0, +4
	36	-5, -2.5, +0, +2.5, +5, +7.5
delta Biolox C-Taper (when used with C-Taper Sleeve - catalog #17-0000E)	28	-2.5, +0, +2.5, +5
	32	-2.5, +0, +2.5, +5
	36	-5, -2.5, +0, +2.5, +5, +7.5
delta Biolox Universal Taper (when used with Universal Taper Sleeve - catalog #6519-T-XXX)	28	-2.5, +0, +4
	32	-2.5, +0, +4
	36	-2.5, +0, +4
	40	-2.5, +0, +4
	44	-2.5, +0, +4

Impact Head onto Calcar Body Trunnion

Select the appropriate corresponding V40 Femoral Head (CoCr, Alumina Ceramic, Biolox delta Ceramic) or sleeve and place it onto the dry trunnion of the femoral stem with a slight twist. Impact the the head with two moderate impactions using the Femoral Head Impactor (Figures 38 and 39).

Reduce Joint and Close

Relocate the Femoral Head into the acetabular cup and check the stability and range of motion. The surgical site is then closed according to the standard procedure for the surgical approach chosen.

Postoperative Care

Postoperative care should progress according to surgeon preference and recommendation.

Figure 37




Figure 38




Figure 39

Optional Step

When selecting a BIOLOX delta Universal Taper Ceramic Femoral Head (6519-1-0xx) for implantation, use of a Universal Adaptor Sleeve is necessary.

Catalog No.	Offset (mm)	Taper
6519-T-025	-2.5	V40
6519-T-100	+0	V40
6519-T-204	+4	V40

After completing the trialing process, intra-operatively assemble the adaptor sleeve to the femoral stem manually. The Universal Adaptor Sleeve must be fully seated on the stem taper before the head is assembled.

Note

In no instance should any attempt be made to pre-assemble the adaptor sleeve inside the BIOLOX delta Universal Ceramic Head.

Intra-operatively assemble the BIOLOX delta Universal Taper Ceramic Head onto the sleeved femoral stem and set with two moderate strikes using the Head Impactor. Care must be taken to avoid excessive impact forces when assembling the Ceramic Head to the sleeved femoral component.

Restoration® Modular

Surgical Protocol

Restoration Modular Calcar Body/Conical Distal Stem Removal

If new components are to be disassembled during surgery (i.e., to readjust version) inspect the proximal body and distal stem closely for damage prior to re-impacting the body onto the distal stem. If the proximal body or distal stem shows damage, do not reuse the components but instead re-implant new, undamaged components.

Note: The Locking Bolt must be removed prior to using stem removal instruments (Figure 40).

Calcar Body Removal

The Body/Stem Separator is made up of three parts: Jackscrew, Shaft Puller, and a reverse-thread Distal Collet (Figure 41). Two modular handles are also available for use with the Body/Stem Separator, which assist in counter-rotation when tightening with the T-Handle.

Unthread the Jackscrew completely from the Shaft Puller prior to inserting through the Calcar Body. Ensure that the Distal Collet is fully threaded into the Shaft Puller, keeping in mind that the Collet and Shaft Puller are reverse-threaded. Insert the Shaft Puller/Distal Collet assembly through the Calcar Body until the collet is fully inserted. An audible click will be heard along with a decrease in resistance upon full insertion.

Thread the Jackscrew through the Shaft Puller/Distal Collet by hand until the Jackscrew cannot be advanced further. Insert the modular handle(s) into the upper hub of the Shaft Puller. The handles are spring-loaded and will engage when rotated to the correct position. Assemble the T-Handle to the Jackscrew and turn the T-Handle until the Calcar Body disengages from the distal stem (Figure 42).

Note: In order to remove the body from the Shaft Puller assembly, remember that the Distal Collet is a *REVERSE THREAD*, and must be completely removed from the assembly to release the body.

Figure 41



Figure 40






Figure 42



Restoration® Modular

Surgical Protocol



Restoration Modular Calcar Body/Conical Distal Stem Removal (continued)

Distal Stem Removal

Assemble the Distal Stem Removal Adapter to the McReynolds Driver-Extractor. Thread the distal stem removal assembly into the insertion feature of the Conical Distal Stem (Figure 43). Use the slap hammer to remove the Conical Distal Stem from the canal.

Removal of the Restoration Modular Calcar Body/ Conical Distal Stem Assembly

The Distal Stem Removal Adapter/McReynolds Driver-Extractor assembly may be threaded through the Calcar Body into the distal stem to remove the entire stem assembly. Use the slap hammer to remove the stem assembly from the canal (Figures 44 and 45).

Figure 43





Figure 44

Figure 45



Restoration® Modular

Surgical Protocol

Calcar Bodies



Calcar Bodies	
CATALOG NO.	SIZE
6276-4-119	19mm +10mm
6276-4-219	19mm +20mm
6276-4-319	19mm +30mm
6276-4-121	21mm +10mm
6276-4-221	21mm +20mm
6276-4-321	21mm +30mm
6276-4-123	23mm +10mm
6276-4-223	23mm +20mm
6276-4-323	23mm +30mm
6276-4-125	25mm +10mm
6276-4-225	25mm +20mm
6276-4-325	25mm +30mm
6276-4-127	27mm +10mm
6276-4-227	27mm +20mm
6276-4-327	27mm +30mm
6276-4-129	29mm +10mm
6276-4-229	29mm +20mm
6276-4-329	29mm +30mm
6276-4-131	31mm +10mm
6276-4-231	31mm +20mm
6276-4-331	31mm +30mm

Restoration® Modular

Surgical Protocol

Conical Distal Stems



Conical Distal Stems-Straight

CATALOG NO.	SIZE
6276-7-014	14mm x 155mm
6276-7-015	15mm x 155mm
6276-7-016	16mm x 155mm
6276-7-017	17mm x 155mm
6276-7-018	18mm x 155mm
6276-7-019	19mm x 155mm
6276-7-020	20mm x 155mm
6276-7-021	21mm x 155mm
6276-7-022	22mm x 155mm
6276-7-023	23mm x 155mm
6276-7-024	24mm x 155mm
6276-7-025	25mm x 155mm
6276-7-026	26mm x 155mm
6276-7-027	27mm x 155mm
6276-7-028	28mm x 155mm

Conical Distal Stems-Straight

CATALOG NO.	SIZE
6276-7-114	14mm x 195mm
6276-7-115	15mm x 195mm
6276-7-116	16mm x 195mm
6276-7-117	17mm x 195mm
6276-7-118	18mm x 195mm
6276-7-119	19mm x 195mm
6276-7-120	20mm x 195mm
6276-7-121	21mm x 195mm
6276-7-122	22mm x 195mm
6276-7-123	23mm x 195mm
6276-7-124	24mm x 195mm
6276-7-125	25mm x 195mm
6276-7-126	26mm x 195mm
6276-7-127	27mm x 195mm
6276-7-128	28mm x 195mm

Conical Distal Stems-Bowed

CATALOG NO.	SIZE
6276-7-214	14mm x 195mm
6276-7-215	15mm x 195mm
6276-7-216	16mm x 195mm
6276-7-217	17mm x 195mm
6276-7-218	18mm x 195mm
6276-7-219	19mm x 195mm
6276-7-220	20mm x 195mm
6276-7-221	21mm x 195mm
6276-7-222	22mm x 195mm
6276-7-223	23mm x 195mm
6276-7-224	24mm x 195mm
6276-7-225	25mm x 195mm
6276-7-226	26mm x 195mm
6276-7-227	27mm x 195mm
6276-7-228	28mm x 195mm


Conical Distal Stems-Bowed

CATALOG NO.	SIZE
6276-7-314	14mm x 235mm
6276-7-315	15mm x 235mm
6276-7-316	16mm x 235mm
6276-7-317	17mm x 235mm
6276-7-318	18mm x 235mm
6276-7-319	19mm x 235mm
6276-7-320	20mm x 235mm
6276-7-321	21mm x 235mm
6276-7-322	22mm x 235mm
6276-7-323	23mm x 235mm
6276-7-324	24mm x 235mm
6276-7-325	25mm x 235mm
6276-7-326	26mm x 235mm
6276-7-327	27mm x 235mm
6276-7-328	28mm x 235mm

Restoration® Modular

Surgical Protocol

Restoration® Modular Instrument System Starter Tray #6278-9-900	
6278-1-150	Resection Guide
6278-5-200	Starter Awl
6278-5-250	Box Chisel
6278-9-090	Large T-Handle
1101-2100	Small T-Handle
1120-1000	Mallet
6278-5-300	Clear Out Reamer




Restoration® Modular

Surgical Protocol


Conical Distal Reamer Tray #1 13mm - 20mm #6278-9-910

6278-8-213	13mm Conical Distal Reamer
6278-8-214	14mm Conical Distal Reamer
6278-8-215	15mm Conical Distal Reamer
6278-8-216	16mm Conical Distal Reamer
6278-8-217	17mm Conical Distal Reamer
6278-8-218	18mm Conical Distal Reamer
6278-8-219	19mm Conical Distal Reamer
6278-8-220	20mm Conical Distal Reamer



Conical Distal Reamer Tray #2 21mm - 28mm #6278-9-911


6278-8-221	21mm Conical Distal Reamer
6278-8-222	22mm Conical Distal Reamer
6278-8-223	23mm Conical Distal Reamer
6278-8-224	24mm Conical Distal Reamer
6278-8-225	25mm Conical Distal Reamer
6278-8-226	26mm Conical Distal Reamer
6278-8-227	27mm Conical Distal Reamer
6278-8-228	28mm Conical Distal Reamer



Restoration® Modular

Surgical Protocol

Proximal Cone Reamer Tray #6278-9-942	
6278-1-519	19mm Proximal Cone Reamer
6278-1-521	21mm Proximal Cone Reamer
6278-1-523	23mm Proximal Cone Reamer
6278-1-525	25mm Proximal Cone Reamer
6278-1-527	27mm Proximal Cone Reamer
6278-1-529	29mm Proximal Cone Reamer
6278-1-531	31mm Proximal Cone Reamer
6278-9-500	Proximal Cone Reamer Post




Restoration® Modular

Surgical Protocol

Calcar Instrument Tray #6278-9-962

6278-4-500	Calcar Body Resection Guide
6278-4-550	Calcar Body Drill Guide
6278-4-400	Calcar Post Resection Cutting Block
6278-4-401	Calcar Post Resection Ring - R
6278-4-402	Calcar Resection Guide Post
6278-4-403	Calcar Post Resection Ring - L
6278-4-450	Lateral Drill Guide
5085-1-227	2.2mm x 150mm Drill (2 per tray)
5800-5-125	3.2mm x 125mm Drill (2 per tray)




Calcar Body Trial Tray #1 19mm - 25mm #6278-9-960

6278-0-119	19mm +10mm Calcar Body Trial
6278-0-219	19mm +20mm Calcar Body Trial
6278-0-319	19mm +30mm Calcar Body Trial
6278-0-121	21mm +10mm Calcar Body Trial
6278-0-221	21mm +20mm Calcar Body Trial
6278-0-321	21mm +30mm Calcar Body Trial
6278-0-123	23mm +10mm Calcar Body Trial
6278-0-223	23mm +20mm Calcar Body Trial
6278-0-323	23mm +30mm Calcar Body Trial
6278-0-125	25mm +10mm Calcar Body Trial
6278-0-225	25mm +20mm Calcar Body Trial
6278-0-325	25mm +30mm Calcar Body Trial



Calcar Body Trial Tray #2 27mm - 31mm #6278-9-961

6278-0-127	27mm +10mm Calcar Body Trial
6278-0-227	27mm +20mm Calcar Body Trial
6278-0-327	27mm +30mm Calcar Body Trial
6278-0-129	29mm +10mm Calcar Body Trial
6278-0-229	29mm +20mm Calcar Body Trial
6278-0-329	29mm +30mm Calcar Body Trial
6278-0-131	31mm +10mm Calcar Body Trial
6278-0-231	31mm +20mm Calcar Body Trial
6278-0-331	31mm +30mm Calcar Body Trial





Restoration® Modular

Surgical Protocol

Finishing Instrument Tray #1 #6278-9-970

6278-1-100	Version Control Stem Inserter
6278-1-110	Stem Inserter Proximal Body Impactor
6278-9-070	Body/Stem Separator
8000-0000	Body/Stem Separator Handles (2 per tray)
6278-9-080	McReynolds Proximal Body Adapter
6266-0-140	Head Impactor
6260-4-070	Proximal Body Steady Handle
6260-4-080	Torque Wrench Adapter
6260-4-090	McReynolds Distal Stem Adapter
6278-5-100	5mm Hex Locking Bolt Driver
6278-5-120	8mm Hex Locking Bolt Driver
6060-2-640	Torque Wrench
6264-8-122R	22mm +0 Head Trial
6264-8-222R	22mm +3 Head Trial
6264-8-322R	22mm +8 Head Trial
6264-8-026R	26mm -3 Head Trial
6264-8-126R	26mm +0 Head Trial
6264-8-226R	26mm +4 Head Trial
6264-8-326R	26mm +8 Head Trial
6264-8-426R	26mm +12 Head Trial
6264-8-028R	28mm -4 Head Trial
6264-8-128R	28mm +0 Head Trial
6264-8-228R	28mm +4 Head Trial
6264-8-328R	28mm +8 Head Trial
6264-8-428R	28mm +12 Head Trial
6264-8-032R	32mm -4 Head Trial
6264-8-132R	32mm +0 Head Trial
6264-8-232R	32mm +4 Head Trial
6264-8-332R	32mm +8 Head Trial
6264-8-432R	32mm +12 Head Trial
6264-8-036R	36mm -5 Head Trial
6264-8-136R	36mm +0 Head Trial
6264-8-236R	36mm +5 Head Trial
6264-8-336R	36mm +10 Head Trial



Finishing Instrument Tray #2 #6278-9-971

6869-1-000	Shaft, McReynolds Extractor Assembly
6869-2-000	Driving Handle, McReynolds Extractor Assembly
6869-3-000	Sliding Hammer, McReynolds Extractor Assembly
6278-1-200D	Distal Stem Impactor
6278-1-350	Proximal Body Impactor



Reconstructive

Hips
Knees
Trauma & Extremities
Foot & Ankle
Joint Preservation
Orthobiologics & Biosurgery

MedSurg

Power Tools & Surgical Accessories
Computer Assisted Surgery
Endoscopic Surgical Solutions
Integrated Communications
Beds, Stretchers & EMS
Reprocessing & Remanufacturing

Neurotechnology & Spine

Craniomaxillofacial
Interventional Spine
Neurosurgical, Spine & ENT
Neurovascular
Spinal Implants

**325 Corporate Drive
Mahwah, NJ 07430
t: 201 831 5000**

www.stryker.com

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. The products depicted are CE marked according to the Medical Device Directive 93/42/EEC. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Dall-Miles, PureFix, Restoration, Stryker, Stryker Orthopaedics, V40. All other trademarks are trademarks of their respective owners or holders.

RMOD-SP-10
2/15

Copyright © Stryker 2015
Printed in USA.