

Triathlon® PKR

Partial Knee Resurfacing Surgical Protocol



Triathlon PKR

Surgical Protocol

Table of Contents

Acknowledgments	4
Introduction	4
Procedural Summary/Highlights	5
Tibial Preparation	8
Excision of Osteophytes	8
Assembling the Ankle Clamp	9
Attaching the Ankle Clamp	10
Flexion/Extension Alignment (Tibial Slope)	11
Varus/Valgus Adjustment	11
Tibial Resection Level	12
Sagittal Resection and Rotational Alignment	13
Transverse Tibial Resection	13
Before Removing the Tibial Resection Guide	13
Femoral Preparation	15
Distal Resection	15
Alignment Check	16
Posterior and Chamfer Resections	17
Sizing and Trialing	19
Tibial Sizing	19
Femoral Trialing	20
Trial Reduction	21
Femoral Peg Preparation	21
Tibial Peg Preparation	22
Cementation	23
Tibial Cementation	23
Femoral Cementation	23
Tibial Insert	24
Alternate Femoral Preparation	25



Triathlon PKR

Surgical Protocol

Acknowledgments

Special thanks to all those surgeons involved in the design and development of the Triathlon PKR System.

Introduction

Triathlon Partial Knee Resurfacing is a fixed bearing unicompartmental knee that blends the ease of a measured resection system with the control of a ligament tension restoring gap balancing system. The single radius design of the Triathlon PKR implant is designed to produce more stable, natural knee motion by coupling the ability to restore ligament tension and balance flexion and extension gaps.

- ▶ The femoral implant has a single radius from 10° of flexion through 110° of flexion and is equally 7mm thin posteriorly and distally. This design allows for gap balancing in only a few simple steps.
- ▶ The Triathlon PKR balances flexion and extension gaps by varying the distal resection. In the Triathlon PKR System the posterior resection remains a fixed 7mm, while the distal resection is varied to accommodate degenerative changes of the condyle.



Procedural Summary/Highlights

1 Tibial Resection



2 Assess Flexion Gap



3 Assess Extension Gap



4 Verify Alignment



5 Distal Resection



6 Posterior/Chamfer Resection



7 Trialing



8 Final Preparation



9 Cementation



Triathlon PKR

Surgical Protocol

Indications

The Triathlon PKR components are for use in unicompartmental knee arthroplasty as a result of:

- Moderately disabling joint disease of the knee resulting from painful osteo- or post traumatic arthritis.
- Revision of previous unsuccessful surgical procedures, either involving, or not involving, previous use of a unicompartmental knee prosthesis.
- As an alternative to tibial osteotomy in patients with unicompartmental osteoarthritis.
- Where bone stock is of poor quality or inadequate for other reconstructive techniques as indicated by deficiencies of the femoral condyle/tibial plateau.

These components are intended for implantation with bone cement.

Contraindications include:

- Patient has an active or suspected latent infection in or about the knee joint.
- Patient has a known sensitivity to device materials.
- Patient's bone stock is compromised by disease and/or infection, or prior implantation which cannot provide adequate support and/or fixation cannot be provided to the prosthesis.
- Patients with inflammatory arthritis.
- Patients with major deformity affecting the mechanical axis of the knee or neuromuscular disorders compromising motor control and/or stability.
- Any mental neuromuscular disorder, which would create an unacceptable risk of prosthesis instability, prosthesis fixation failure or complications in post-operative care.
- Skeletal immaturity.
- Ligamentous instability such that the post-operative stability afforded by the unicompartmental knee prosthesis would be compromised such as multidirectional/ACL instability.
- Untreated damage to the contralateral compartment or the ipsilateral knee not being replaced by a prosthesis.
- Untreated deterioration or destruction of the patello-femoral joint.
- Severe deformity and/or recurrent subluxation of the knee joint.
- Obesity. An overweight or obese patient can produce loads on the prosthesis, which can lead to failure of fixation of the device or failure of the device itself.
- Severe tibial bone loss/deformity (over 15 degrees varus).

* See package insert for warnings, precautions, adverse effects and other essential product information.

Surgical Procedure

Triathlon PKR

Surgical Protocol



Figure 1

Tibial Preparation

Triathlon PKR System utilizes a “**tibia first**” approach to the unicompartmental arthroplasty. Accurate alignment of instrumentation to anatomical references during the tibial resection step of this procedure will aid in establishing tibial/femoral component alignment. Since all femoral resections are instrumented by referencing the tibial resection, care should be taken during the preparation of the tibial bone.




Figure 2

Excision of Osteophytes

- ▶ Before beginning bone preparation of the affected compartment, there are several areas of the joint that should be inspected for osteophytes. All osteophytes must be removed before proceeding to the rest of the procedure. These areas include but are not limited to the:
 - ▶ Margin of the femoral condyle.




Figure 3

- ▶ Margins of the intercondylar notch.



Figure 4

- ▶ Beneath the medial collateral ligament.
- ▶ Additional preparation of the affected compartment may also include the removal of the meniscus.



Assembling the Ankle Clamp

Figure 5

Triathlon PKR

Surgical Protocol



Figure 6

Attaching the Ankle Clamp

- ▶ The fixation arm of the Proximal Rod is fully extended to reach the tibial eminence. Pin the head of the Proximal Rod into an exposed area of the anterior proximal tibia, generally located anterior to the ACL insertion point.




Figure 7

Note: Fully extending the fixation arm prior to attaching it to the Tibial eminence will help allow for clearance of the Tibial Resection Guide. Once the headless pin has been placed, the Tibial resection guide may be advanced toward the tibia.

- ▶ Move the leg to 70° of flexion (flexion past 70° may result in tightening of the quadriceps muscle and make soft tissue retraction difficult).
- ▶ Attach the Ankle Clamp assembly to the ankle and position the distal tube of the assembly along the long access of the tibia.
- ▶ Rotate the entire assembly to ensure that the base of the assembly is aligned with the center of the ankle. The center of the ankle is generally in line with the second metatarsal or medial 1/3rd of the tibial tubercle.



Figure 8



Figure 9

Flexion/Extension Alignment (Tibial Slope)

- ▶ Flexion/Extension alignment is correct when the long axis of the Ankle Clamp/Proximal Rod Assembly parallels the mid-coronal plane of the tibia. This scenario will provide 0° of slope.

Note: The slope of the Tibial Resection Guide is 0° degrees. There is no slope built into the tibial components of the Triathlon PKR system.

- ▶ Additional slope can be introduced to the guide by pulling the distal end of the Ankle Clamp Assembly away from the ankle.



Figure 10

Varus/Valgus Adjustment

- ▶ Varus/Valgus Alignment can be changed by loosening the adjustment knob on the Ankle Clamp Assembly and by sliding the assembly medially/laterally.
- ▶ Once proper alignment is achieved, re-tighten all adjustment knobs on the Ankle Clamp Assembly.

Triathlon PKR

Surgical Protocol



Figure 11

Tibial Resection Level

- ▶ Attach the Tibial Stylus to the Tibial Resection Guide. With the “4” end referencing the lowest level of the tibial plateau. 4mm of bone will be resected below the tip of the stylus. Alternatively, if the “2” end of the Tibial Stylus is used, the amount of bone resected will be 2mm below the tip of the stylus.



Figure 12

- ▶ The height of the Tibial Stylus and Tibial Resection Guide can be adjusted using the wheel on the Tibial Resection Guide.
- ▶ For coarse adjustment, press the wheel between thumb and forefinger and slide the assembly.
- ▶ For fine adjustment, turn the wheel counter-clockwise to decrease the tibial resection level or turn clockwise to increase the tibial resection level.
- ▶ If additional stability of the guide is required, a pin may be inserted through the fixation hole below the tibial resection guide.




Figure 13



Figure 14

Sagittal Resection and Rotational Alignment

- ▶ The Sagittal resection (vertical tibial resection), is made with a reciprocating saw flush against the tibial spine. In the medial compartment, the rotation of the resection is estimated by rotating the saw blade such that it points toward the center of the femoral head. The sagittal resection should be made just medial to the origin of the ACL avoiding damage to its fibers. The sagittal resection is complete when the saw comes to a stop on the upper surface of the Tibial Resection Guide.

Note: Before making the transverse resection, a retractor should be inserted to protect the medial collateral ligament.

Transverse Tibial Resection

- ▶ This resection should be continued all of the way to the back of the joint to ensure a consistent flat and smooth tibial surface. When the plateau is loose it is levered up with a broad osteotome and removed. Soft tissue attachments posteriorly may need to be cut with a knife. The posterior horn of the medial meniscus can now be removed, if present. Care should be taken to avoid saw blade excursion to extend beyond the intersection of the sagittal resection and under the tibial eminence.



Figure 15

Before Removing the Tibial Resection Guide

- ▶ Ankle Clamp Assembly, the flexion gap must be checked to ensure that enough tibial bone has been resected to accommodate an 8mm tibial component (the thinnest tibial component offered with this system).
- ▶ Remove retractors from the joint space (if they are left in, they have the effect of tightening the soft tissues which artificially diminishes the gap) and insert the 8mm spacer block. If it cannot be inserted, or feels too tight, then more bone needs to be excised from the tibia.
- ▶ To make a 2mm skim cut from the tibia, the stylus is re-assembled to the tibial resection guide. The wheel on the Tibial Resection Guide is then adjusted until the 2mm stylus tip touches the prepared surface of the tibial plateau. The tibial resection process is then repeated. An additional 2mm of tibial bone will be removed.

Triathlon PKR

Surgical Protocol

Tibial Preparation

Femoral Preparation




Figure 16

- ▶ Assembly of Alignment Guide.
- ▶ Assembly of Distal Resection Guide.
- ▶ Assembly of 2:1 Cutting Block.

Locked



Figure 17

Unlocked



Figure 18



Figure 19




Figure 20



Figure 21

Femoral Preparation

Distal Resection

- ▶ The formula for determining which distal resection guide to use is as follows:

Flexion Gap – Extension Gap = Distal Resection Guide

- ▶ Determine the appropriate Distal Resection Guide

Remove all retractors from the joint.

- 1) Assess flexion gap using the Spacer Block Handles.
- 2) Assess extension gap using the Spacer Block Handles.
- 3) Flexion Gap (mm) – Extension Gap (mm) = Distal Resection Guide.

Example A

8mm – 10mm = -2 or “Minus 2” Distal Resection Guide

Example B


9mm – 10mm = -1 or “Minus 1” Distal Resection Guide

Note: See page 25 for Gap Balancing reference chart.

- ▶ Attach the appropriate Distal Resection Guide to the Spacer Block Handle equal to the assessed extension gap. The open side of the resection block should be facing the outside of the joint when the assembly is inserted into the joint.
- ▶ With the leg in extension, insert the Spacer Block Handle into the extension gap with Distal Resection Guide attached.
- ▶ Ensure that the Spacer Block Handle is inserted far enough into the joint so that it is completely under the distal femoral bone.

Triathlon PKR

Surgical Protocol



Alignment Check

Note: From the lateral view, confirm a neutral distal resection. If slope has been introduced into the tibial resection, the spacer block will press against the femur in an extended position. To achieve neutral flexion of the distal resection, flex the tibia until the alignment rod is parallel to the long axis of the femur in the sagittal plane. A rolled towel or object placed under the knee may aid in this process. Pin the Distal Resection Guide into place once proper alignment is achieved (neutral distal resection).

Figure 22






Figure 23

- ▶ Before performing the distal resection, retraction of the MCL is recommended to prevent saw excursion damage.
- ▶ Perform the distal femoral resection. Do not extend the saw blade posteriorly past the distal femur to avoid damaging the posterior popliteal region. If desired, the distal resection can be started in extension and finished in flexion. The Spacer Block must be removed from the joint before the knee is placed in flexion.

Posterior and Chamfer Resections

- ▶ Attach the 2-in-1 Cutting Block to the Spacer Block Handle that fit into the determined (assessed) flexion gap.

Note: The system contains a left medial/right lateral and right medial/left lateral 2-in-1 Cutting Block. The cutting blocks are marked for usage on the appropriate condyle.



Figure 24

- ▶ Flex the Knee to 90° Degrees.
- ▶ Insert the Spacer Block Handle with the 2-in-1 Cutting Block attached into the flexion gap. The Spacer Block Handle must be fully inserted so that it extends beyond the most posterior aspect of the condyle. This will help insure the proper posterior resection level (7mm).
- ▶ Adjust the flexion and tibial rotation of the knee so that the Spacer Block Handle rests flat on the tibial resection and back side of the **2-in-1 Cutting Block rests flat on the resected distal femur.**




Figure 25

- ▶ Proper placement of the 2-in-1 Cutting Block is achieved when the Spacer Block Handle contacts both the posterior condyle and the tibial resection and the **2-in-1 Cutting Block rests flat on the resected distal femur.** The 2-in-1 Cutting Block should be positioned so the anterior most anterior hole in the guide is centered on the distal resection.
- ▶ Pin the 2-in-1 Cutting Block into place.
- ▶ Unlock and remove Spacer Block.
- ▶ Perform posterior resection.

Triathlon PKR

Surgical Protocol

Femoral
Preparation



Figure 26

- ▶ Perform chamfer resection.

Sizing and
Trialing




Figure 27

- ▶ Check for and remove posterior condylar osteophytes if necessary.



Figure 28

Sizing and Trialing

Tibial Sizing

- Using the Tibial Sizer, determine the best fit tibial component size. The Tibial Sizer should be placed flat on the transverse tibial resection and against the wall of the sagittal resection. The lip on the sizer should be hooked over the posterior edge of the tibia and pulled forward toward the surgeon. A secondary check of the tibial size can be made by checking the medial-lateral width of the tibial resection using the opposite end of the sizer at the widest cross-section of the tibial resection. If the A/P and M/L measurements do not both match the same size tibial component, the surgeon may consider increasing the tibial sagittal resection to improve tibial component fit.



Figure 29

M/L Check



Figure 30

Triathlon PKR

Surgical Protocol

Incorrect




Figure 31

Correct




Figure 32

Femoral Trialing

- ▶ Bring the knee into approximately 90 degrees of flexion.
- ▶ Select the Femoral Trial that most closely approximates the dimensions of the prepared femur. Insert the Femoral Trial into the incision and assemble onto the femoral condyle. Position the trial so that it rests on the distal and posterior surfaces.
- ▶ A properly sized femoral component should have 2mm-3mm of exposed bone above the anterior edge of the trial.

Note: The flexibility of this system allows for one up, one down, femoral and tibial size interchangeability.



Figure 33

- ▶ Repeat with additional trials until the proper size is selected. If the femur appears to be between two sizes, choose the smaller size. This will help to ensure proper patella tracking.
- ▶ Once the appropriate femoral component is chosen, insert the Tibial Trial beneath the femoral component.

Trial Reduction

- ▶ Trial reduction allows the surgeon to check for accurate tibio-femoral tracking.
- ▶ A final check is made, with valgus and varus stress, to ensure that desired degree of varus/valgus laxity has been achieved.



Figure 34

Femoral Peg Preparation

- ▶ Once the appropriate size and location of the femoral component is determined, the Femoral Trial can be pinned anteriorly for additional stability while preparing for femoral pegs.




Figure 35

- ▶ The Femoral Trial Drill Guide is attached to the Modular Handle and assembled into the holes on the articulation surface of the Femoral Trial. Using the appropriately sized Peg Drill (Small Drill for sizes 1-2, Large Drill for sizes 3-6), both holes are drilled. The Peg Drill is advanced until the step on the drill makes contact with the front face of the Femoral Trial Drill Guide.

Note: To prevent possible impingement of the drill, fully engage the Peg Drill in the Femoral Trial Drill Guide before beginning to drill for the femoral peg holes.

Triathlon PKR

Surgical Protocol



Figure 36

- ▶ The Femoral Trial Drill Guide is removed from the trial and the Peg Trial is pushed into the posterior hole of the Femoral Trial. The knee is taken through a final range of motion to confirm component placement and size.
- ▶ The Femoral Trial Peg is removed with the Modular handle. Both the Femoral Trial and Tibial Trial are removed.




Figure 37

Tibial Peg Preparation

- ▶ Tibial preparation is completed by inserting the appropriately sized Tibial Template into the joint. The Tibial Template should be flush with the wall of the sagittal resection and the anterior contour of the template even with the anterior cortical rim of the tibia.
- ▶ The template is held into place with spikes on the bottom of the template. A headless pin may be used for additional stabilization.
- ▶ Using the appropriately sized Peg Drill (Small Drill for sizes 1-2, Large Drill for sizes 3-6), both holes are drilled. The Peg Drill is advanced until the step on the drill makes contact with the face of the Tibial Template.



Figure 38



Figure 39



Figure 40

Cementation

Tibial Cementation

- ▶ Cement should be applied to both bone and implant surfaces to insure the best interdigitation.
- ▶ Press cement into the bone of the tibia, including the wall of the tibial eminence.
- ▶ Apply cement to the bottom and side wall of the component. Minimize the amount of cement used posteriorly on the implant.
- ▶ Impact tibial component into place. Remove excess cement from edge of tibia. Care should be taken to minimize excess cement in posterior aspect of the tibia.

Note: The interior area and rim of the tibial tray must be free of cement and debris to ensure a secure lock between the tibial insert and tibial tray.

Femoral Cementation

- ▶ Cement should be applied to both bone and implant surfaces to ensure the best interdigitation.
- ▶ Press cement into the bone of the femur.
- ▶ Apply cement to the underside of the femoral component. Minimize the amount of cement used posteriorly on the implant.
- ▶ Impact Femoral component into place. Remove excess cement from edge of component. Care should be taken to minimize excess cement in posterior aspect of the femur.

Triathlon PKR

Surgical Protocol



Figure 41

Tibial Insert

- ▶ Insert tibial insert onto the top of the tibial tray.
- ▶ Push the insert posteriorly into the tray until it engages on the posterior lip of the tibial tray.
- ▶ Using the tibial impactor, tap the front edge of the insert down until it seats in the tibial tray.

Note: The tibial insert has been designed so that there is a small space between the tibial insert and outer rim of the tibial tray. When the insert is properly seated, this space is of uniform thickness around the entire perimeter of outer rim of the tibial tray.

- ▶ The leg is brought out into extension for cement pressurization and again flexed for final cement removal.

Note: Avoid hyper-extension of the leg during the curing of the cement as this may cause posterior lift-off of the tibial component. A thin spacer device may be inserted into the joint before the leg is brought into extension to serve as additional cement pressurization during cement curing.




Figure 42

Gap Balancing

*Flexion minus Extension =
Distal Resection Guide*

		Extension Space (mm)							
		8	9	10	11	12	13	14	15
Flexion Space (mm)	8	Zero	Minus 1	Minus 2	Minus 3				
	9	Plus 1	Zero	Minus 1	Minus 2	Minus 3			
	10	Plus 2	Plus 1	Zero	Minus 1	Minus 2	Minus 3		
	11		Plus 2	Plus 1	Zero	Minus 1	Minus 2	Minus 3	
	12			Plus 2	Plus 1	Zero	Minus 1	Minus 2	Minus 3

Cementation

Triathlon PKR

Surgical Protocol

Catalog #	Description	Quantity in Kit
Triathlon PKR Tibial & Femoral Prep Tray (Lower) Part Numbers		
6541-4-602	Universal Alignment Rods	1
6633-7-605	Headless Pin Puller	1
6541-4-401	Navigation Tracker Adapter	Optional
6003-060-000	Resection Plane Probe	Optional
5650-2-600	Tibial Sizer	1
5650-2-601	Tibial Template #1 LM/RL	1
5650-2-602	Tibial Template #2 LM/RL	1
5650-2-603	Tibial Template #3 LM/RL	1
5650-2-604	Tibial Template #4 LM/RL	1
5650-2-605	Tibial Template #5 LM/RL	1
5650-2-606	Tibial Template #6 LM/RL	1
5650-2-611	Tibial Template #1 RM/LL	1
5650-2-612	Tibial Template #2 RM/LL	1
5650-2-613	Tibial Template #3 RM/LL	1
5650-2-614	Tibial Template #4 RM/LL	1
5650-2-615	Tibial Template #5 RM/LL	1
5650-2-616	Tibial Template #6 RM/LL	1
5650-8-101	PKR Tibial & Femoral Prep Tray (Lower)	1
		Total Quantity 16

Catalog #	Description	Quantity in Kit
Triathlon PKR Tibial & Femoral Prep Tray (Upper) Part Numbers		
7650-1038A	3.5" Fluted Headless Pins	1
7650-1039A	2.5" Fluted Headless Pins	1
7650-1035	Pin Driver	1
5650-0809	Spacer Block Handle 8/9	1
5650-1011	Spacer Block Handle 10/11	1
5650-1213	Spacer Block Handle 12/13	1
5650-1415	Spacer Block Handle 14/15	1
5650-4-806	Spacer Block Alignment Handle	1
5650-4-536A	Peg Drill 3-6	1
5650-4-512A	Peg Drill 1-2	1
8000-1040	Ankle Clamp Assembly	1
5650-2-617	Proximal Rod	1
5650-1-724	Distal Resection Guide, Minus 3	1
5650-1-725	Distal Resection Guide, Minus 2	1
5650-1-726	Distal Resection Guide, Minus 1	1
5650-1-727	Distal Resection Guide, Zero	1
5650-1-728	Distal Resection Guide, Plus 1	1
5650-1-729	Distal Resection Guide, Plus 2	1
5650-1-700	2-in-1 Cutting Block, LM/RL	1
5650-1-701	2-in-1 Cutting Block, RM/LL	1
5650-2-424	Tibial Stylus	1
5650-2-700	Tibial Resection Guide LM/RL Captured	1
5650-2-701	Tibial Resection Guide RM/LL Captured	1
5650-2-702	Tibial Resection Guide LM/RL Uncaptured	1
5650-2-703	Tibial Resection Guide RM/LL Uncaptured	1
5650-8-001	PKR Tibial & Femoral Prep Tray (Upper)	1
Total Quantity		26

Triathlon PKR

Surgical Protocol

Catalog #	Description	Quantity in Kit
Triathlon PKR Tibial & Femoral Trials Tray (Upper) Part Numbers		
6541-4-710	Posterior Osteophyte Chisel	1
6541-4-810	Impaction Handle	1
6541-4-808	Modular Handle	1
5650-4-805	Tibial Trial Inserter/Extractor	1
5650-4-811	Femoral Impactor	1
5650-4-812	Tibial Impactor	1
5650-1-712	Femoral Trial Drill Guide LM/RL	1
5650-1-713	Femoral Trial Drill Guide RM/LL	1
5630-T-108	Tibial Trial #1 - 8mm	1
5630-T-109	Tibial Trial #1 - 9mm	1
5630-T-110	Tibial Trial #1 - 10mm	1
5630-T-112	Tibial Trial #1 - 12mm	1
5630-T-208	Tibial Trial #2 - 8mm	1
5630-T-209	Tibial Trial #2 - 9mm	1
5630-T-210	Tibial Trial #2 - 10mm	1
5630-T-212	Tibial Trial #2 - 12mm	1
5630-T-308	Tibial Trial #3 - 8mm	1
5630-T-309	Tibial Trial #3 - 9mm	1
5630-T-310	Tibial Trial #3 - 10mm	1
5630-T-312	Tibial Trial #3 - 12mm	1
5630-T-408	Tibial Trial #4 - 8mm	1
5630-T-409	Tibial Trial #4 - 9mm	1
5630-T-410	Tibial Trial #4 - 10mm	1
5630-T-412	Tibial Trial #4 - 12mm	1
5630-T-508	Tibial Trial #5 - 8mm	1
5630-T-509	Tibial Trial #5 - 9mm	1
5630-T-510	Tibial Trial #5 - 10mm	1
5630-T-512	Tibial Trial #5 - 12mm	1
5630-T-608	Tibial Trial #6 - 8mm	1
5630-T-609	Tibial Trial #6 - 9mm	1
5630-T-610	Tibial Trial #6 - 10mm	1
5630-T-612	Tibial Trial #6 - 12mm	1

Catalog #	Description	Quantity in Kit
Triathlon PKR Tibial & Femoral Trials Tray (Upper) Part Numbers		
5610-T-101	Drill Thru Femoral Trial #1 LM/RL	1
5610-T-102	Drill Thru Femoral Trial #1 RM/LL	1
5610-T-201	Drill Thru Femoral Trial #2 LM/RL	1
5610-T-202	Drill Thru Femoral Trial #2 RM/LL	1
5610-T-301	Drill Thru Femoral Trial #3 LM/RL	1
5610-T-302	Drill Thru Femoral Trial #3 RM/LL	1
5610-T-401	Drill Thru Femoral Trial #4 LM/RL	1
5610-T-402	Drill Thru Femoral Trial #4 RM/LL	1
5610-T-501	Drill Thru Femoral Trial #5 LM/RL	1
5610-T-502	Drill Thru Femoral Trial #5 RM/LL	1
5610-T-601	Drill Thru Femoral Trial #6 LM/RL	1
5610-T-602	Drill Thru Femoral Trial #6 RM/LL	1
5650-4-800	Joint Tension Feeler Gage	1
5610-T-012	Peg Trial 1-2	1
5610-T-036	Peg Trial 3-6	1
6541-4-808	Modular Handle	1
6541-4-810	Impaction Handle	1
5650-8-002	PKR Tibial & Femoral Trials Tray (Upper)	1
		Total Quantity 50
6541-9-000	Triathlon Case (Only Two Needed)	2

Triathlon PKR

Surgical Protocol

Catalog #	Description
Triathlon PKR Insert Part Numbers	
5630-G-108	Triathlon PKR Insert X3 #1 LM/RL - 8mm
5630-G-109	Triathlon PKR Insert X3 #1 LM/RL - 9mm
5630-G-110	Triathlon PKR Insert X3 #1 LM/RL - 10mm
5630-G-112	Triathlon PKR Insert X3 #1 LM/RL - 12mm
5630-G-128	Triathlon PKR Insert X3 #1 RM/LL - 8mm
5630-G-129	Triathlon PKR Insert X3 #1 RM/LL - 9mm
5630-G-120	Triathlon PKR Insert X3 #1 RM/LL - 10mm
5630-G-122	Triathlon PKR Insert X3 #1 RM/LL - 12mm
5630-G-208	Triathlon PKR Insert X3 #2 LM/RL - 8mm
5630-G-209	Triathlon PKR Insert X3 #2 LM/RL - 9mm
5630-G-210	Triathlon PKR Insert X3 #2 LM/RL - 10mm
5630-G-212	Triathlon PKR Insert X3 #2 LM/RL - 12mm
5630-G-228	Triathlon PKR Insert X3 #2 RM/LL - 8mm
5630-G-229	Triathlon PKR Insert X3 #2 RM/LL - 9mm
5630-G-220	Triathlon PKR Insert X3 #2 RM/LL - 10mm
5630-G-222	Triathlon PKR Insert X3 #2 RM/LL - 12mm
5630-G-308	Triathlon PKR Insert X3 #3 LM/RL - 8mm
5630-G-309	Triathlon PKR Insert X3 #3 LM/RL - 9mm
5630-G-310	Triathlon PKR Insert X3 #3 LM/RL - 10mm
5630-G-312	Triathlon PKR Insert X3 #3 LM/RL - 12mm
5630-G-328	Triathlon PKR Insert X3 #3 RM/LL - 8mm
5630-G-329	Triathlon PKR Insert X3 #3 RM/LL - 9mm
5630-G-320	Triathlon PKR Insert X3 #3 RM/LL - 10mm
5630-G-322	Triathlon PKR Insert X3 #3 RM/LL - 12mm

Catalog #	Description
-----------	-------------

Triathlon PKR Insert Part Numbers

5630-G-408	Triathlon PKR Insert X3 #4 LM/RL - 8mm
5630-G-409	Triathlon PKR Insert X3 #4 LM/RL - 9mm
5630-G-410	Triathlon PKR Insert X3 #4 LM/RL - 10mm
5630-G-412	Triathlon PKR Insert X3 #4 LM/RL - 12mm
5630-G-428	Triathlon PKR Insert X3 #4 RM/LL - 8mm
5630-G-429	Triathlon PKR Insert X3 #4 RM/LL - 9mm
5630-G-420	Triathlon PKR Insert X3 #4 RM/LL - 10mm
5630-G-422	Triathlon PKR Insert X3 #4 RM/LL - 12mm
5630-G-508	Triathlon PKR Insert X3 #5 LM/RL - 8mm
5630-G-509	Triathlon PKR Insert X3 #5 LM/RL - 9mm
5630-G-510	Triathlon PKR Insert X3 #5 LM/RL - 10mm
5630-G-512	Triathlon PKR Insert X3 #5 LM/RL - 12mm
5630-G-528	Triathlon PKR Insert X3 #5 RM/LL - 8mm
5630-G-529	Triathlon PKR Insert X3 #5 RM/LL - 9mm
5630-G-520	Triathlon PKR Insert X3 #5 RM/LL - 10mm
5630-G-522	Triathlon PKR Insert X3 #5 RM/LL - 12mm
5630-G-608	Triathlon PKR Insert X3 #6 LM/RL - 8mm
5630-G-609	Triathlon PKR Insert X3 #6 LM/RL - 9mm
5630-G-610	Triathlon PKR Insert X3 #6 LM/RL - 10mm
5630-G-612	Triathlon PKR Insert X3 #6 LM/RL - 12mm
5630-G-628	Triathlon PKR Insert X3 #6 RM/LL - 8mm
5630-G-629	Triathlon PKR Insert X3 #6 RM/LL - 9mm
5630-G-620	Triathlon PKR Insert X3 #6 RM/LL - 10mm
5630-G-622	Triathlon PKR Insert X3 #6 RM/LL - 12mm

Catalog #	Description
-----------	-------------

Triathlon PKR Femur Part Numbers

5610-F-101	Triathlon PKR Femur #1 LM/RL
5610-F-201	Triathlon PKR Femur #2 LM/RL
5610-F-301	Triathlon PKR Femur #3 LM/RL
5610-F-401	Triathlon PKR Femur #4 LM/RL
5610-F-501	Triathlon PKR Femur #5 LM/RL
5610-F-601	Triathlon PKR Femur #6 LM/RL
5610-F-102	Triathlon PKR Femur #1 RM/LL
5610-F-202	Triathlon PKR Femur #2 RM/LL
5610-F-302	Triathlon PKR Femur #3 RM/LL
5610-F-402	Triathlon PKR Femur #4 RM/LL
5610-F-502	Triathlon PKR Femur #5 RM/LL
5610-F-602	Triathlon PKR Femur #6 RM/LL

Triathlon PKR Baseplate Part Numbers

5620-B-101	Triathlon PKR Baseplate #1 LM/RL
5620-B-201	Triathlon PKR Baseplate #2 LM/RL
5620-B-301	Triathlon PKR Baseplate #3 LM/RL
5620-B-401	Triathlon PKR Baseplate #4 LM/RL
5620-B-501	Triathlon PKR Baseplate #5 LM/RL
5620-B-601	Triathlon PKR Baseplate #6 LM/RL
5620-B-102	Triathlon PKR Baseplate #1 RM/LL
5620-B-202	Triathlon PKR Baseplate #2 RM/LL
5620-B-302	Triathlon PKR Baseplate #3 RM/LL
5620-B-402	Triathlon PKR Baseplate #4 RM/LL
5620-B-502	Triathlon PKR Baseplate #5 RM/LL
5620-B-602	Triathlon PKR Baseplate #6 RM/LL

Joint Replacements

Trauma, Extremities & Deformities

Craniomaxillofacial

Spine

Biologics

Surgical Products

Neuro & ENT

Interventional Pain

Navigation

Endoscopy

Communications

Imaging

Patient Handling Equipment

325 Corporate Drive
Mahwah, NJ 07430
t: 201 831 5000

www.stryker.com

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker and Triathlon. All other trademarks are trademarks of their respective owners or holders.

Literature Number: **LTPKR-SP Rev. 4**

Copyright © 2011 Stryker
Printed in USA