

# Metal Augment Trial

## Surgical Protocol Addendum

### Femoral Augment Trial Assembly – Metal Augment Trials

If distal or posterior augments are required in a Posterior Stabilized or Total Stabilized setting, assemble the Metal Augment Trial(s) to the surface of the appropriate size Femoral Trial component corresponding to augmentation resections.

**Note:** Metal Augment Trials are not size specific; blocks are universally sized and offered in 5mm, 10mm and 15mm thicknesses.

Assemble the Metal Augment Trial to the most center position of the augmented surface on the Femoral Trial, as seen in **Figures 1B and 2B**.

Ensure the text on the Metal Augment Trial is facing down (as the magnet is strongest on that side), see **Figures 1B and 2B**.

**Impaction Note: If hand impaction is not sufficient, remove Metal Augment Trials before using a force impaction method and re-assemble to avoid detachment of Metal Augment Trials.**

### TS Trialing



Figure 1A

### PS Trialing



Figure 2A



Figure 1B



Figure 2B



Figure 1C




Figure 2C

# Metal Augment Trial

## Surgical Protocol Addendum

**stryker**<sup>®</sup>

325 Corporate Drive  
Mahwah, NJ 07430  
t: 201.831.5000

[www.stryker.com](http://www.stryker.com)

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Triathlon and Stryker. All other trademarks are trademarks of their respective owners or holders.

LSPK48AD2  
Copyright © 2009 Stryker

MS/GS 12/09  
Printed in USA