### stryker

Trauma

# Gamma3 Plus Targeting System



# **Table of Contents**

| Introduction   | 5 |
|--|---|
| Implant Features of the Trochanteric Gamma3 Nail           | 6 |
| Implant Features of the Long Gamma3 Nail                   | 7 |
| Indications/Contraindications for Trochanteric Gamma3 Nail | 8 |
| Indications/Contraindications for Long Gamma3 Nail         | 9 |

### **Operative Technique**

| Implant Selection   | 10 |
|---|----|
| Patient Positioning   | 11 |
| Fracture Reduction  | 12 |
| Incision  | 13 |
| Entry Point   | 13 |
| Preparation of Medullary Canal for Trochanteric Gamma3 Nail | 14 |
| Preparation of Medullary Canal for Long Gamma3 Nail         | 15 |
| One Step Conical Reamer                                     | 16 |
| Assembly of Targeting Device                                | 17 |
| Nail Insertion  | 19 |
| Optional One Shot Device                                    | 20 |
| Lag Screw Insertion   | 21 |
| Lag Screw Fixation  | 26 |
| Distal Screw Locking of the Trochanteric Gamma3 Nail        | 28 |
| Distal Screw Locking of the Long Gamma3 Nail                | 30 |
| End Cap Insertion   | 32 |
| Nail Extension End Caps                                     | 33 |
| Extraction of the Gamma3 Implants                           | 34 |
| Ordering Information – Implants                             | 36 |
| Ordering Information – Instruments                          | 38 |

### Gamma3 Nail

### **Contributing Surgeons**

#### Christopher T. Born, M.D.

Professor of Orthopaedic Surgery Brown University Providence, RI USA

#### James Maxey, M.D.

University of Illinois College of Medicine Peoria, IL USA

#### Robert Probe, M.D.

Department of Orthopaedic Surgery Texas A&M University Health Science Center Temple, TX USA

#### Gilbert Taglang, M.D.

Association Internationale pour L'Osteosyntheses Dynamique (AIOD) Strasbourg France

We gratefully acknowledge and appreciate the contributions to this technique made by:

#### Kevin Luke, M.D.

Parkview Musculoskeletal Institute. Chicago, IL USA

Daniel Horwitz, M.D.

University of Utah Salt Lake City, UT USA

### Anthony Sorkin, M.D.

Rockford Orthopaedics Associates Rockford, IL USA This publication sets forth detailed recommended procedures for using Stryker devices and instruments.

It offers guidance that you should heed, but, as with any such technical guide, each surgeon must consider the particular needs of each patient and make appropriate adjustments when and as required. A workshop training is required prior to first surgery.

#### Note

All bone screws referenced in this material are not approved for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic or lumbar spine.

# Gamma3 Nail Introduction

#### Introduction

The Gamma3 Locking Nail System is based on more than 15 years of Gamma Nail experience and is the third generation of intramedullary short and long Gamma fixation nails.

This evolution is based on the experience and clinical outcomes from surgeons all over the world.

The Gamma3 Plus Targeting System is designed to facilitate Gamma Nail surgery in patients with excess soft tissue.

Improvements in design and function allow for:

- An implant with a 15.5mm proximal diameter, reducing the amount of bone taken prior to nail insertion.
- A Lag Screw with a 10.5mm diameter and improved cutting flutes while maintaining the superior resistance to cut-out of earlier Lag Screw designs.
- 5mm Distal Locking Screws that are compatible with the T2 IM Nail System.
- Titanium Alloy (Ti-6Al-4V).

Another major advantage of this system is the newly designed instrument platform. The intuitive and ergonomic instrument designs are intended to reduce the learning curve of the OR staff, as well as providing for a less invasive technique, where applicable. Additionally, the Gamma3 instruments share the same design platform as the critically acclaimed T2 and S2 IM Nail Systems.

#### Design Features of the Gamma3 System

The Gamma3 Locking Nails come in three neck-shaft angles (120°, 125° and 130°).

All nails use the same Lag Screw, Set Screw, Distal Locking Screws and End Caps.





# Gamma3 Nail Implant Features

#### **Trochanteric Gamma3 Nail**

The anatomical shape of the nail is universal for all indications involving the treatment of intertrochanteric fractures. The nail is cannulated for controlled insertion and features a conical tip for positioning within the medullary canal.

A range of three different neck-shaft angles is available for Lag Screw entry to accommodate variations in femoral neck anatomy.

A single Distal Locking Screw is provided to stabilize the nail in the medullary canal and to help prevent rotation in complex fractures. The oblong hole allows for static or dynamic locking.

### **Technical Specifications**

- **Material:** Type II anodized Titanium alloy (Ti-6Al-4V)
- Nail length: 180mm
- Nail diameter: Proximal 15.5mm, Distal: 11.0mm
- Proximal Nail neck-shaft angles: 120°, 125°, 130°
- M-L bend for valgus curvature: 4 degrees
- End Caps in lengths of: Standard, +5mm and +10mm
- **Distal oblong hole:** 5mm screws allow up to 5mm of dynamic compression
- Lag Screw diameter: 10.5mm
- Lag Screw lengths: 75–120mm in 5mm increments
- Patented Self-Tapping Lag Screw design for high load absorption and easy insertion

Each nail is supplied sterile, packaged together with a Set Screw in one box.

- Asymmetrical depth profile to allow the Lag Screw to slide in the lateral direction only for dynamic bone compression at the fracture site to enhance fracture healing
- A patented self retaining Set Screw to help prevent rotation and to simultaneously allow sliding of the Lag Screw laterally
- Distal Locking Screw diameter: 5mm
- Fully threaded Distal Locking Screw lengths: ranging from 25-120mm in 2.5mm or 5mm increments
- Self-tapping screw tip with optimized short cutting flutes
- **Optimized diameter** under the head helps to prevent microfractures during insertion
- Low profile head





### **Distal Locking Options**

- Locking in the distal portion of the oblong hole creates a dynamic locking mechanism (Figure 1).
- Locking in the proximal part of the oblong hole allows static locking of the nail (Figure 2).



Figure 1 Dynamic Locking

9

Figure 2 Static Locking

# Gamma3 Nail Implant Features

#### Long Gamma3 Nail

This nail incorporates several important mechanical design features, including an unslotted and cannulated design for Guide-Wire controlled insertion.

The three neck-shaft angles: 120°, 125° and 130° accommodate variations in femoral neck anatomy. The Long Nail offers the opportunity to use two distal Locking Screws that are inserted through the distal end of the nail to control rotation and telescoping. As shown below, the nail offers the possibility for either static, dynamic or secondary dynamic distal locking, depending on the fracture pattern.

### **Technical Specifications**

- **Material:** Type II anodized Titanium alloy (Ti-6Al-4V)
- Nail length: 280mm to 480mm, in 20mm increments, in left and right (shorter or longer nails are available on request)
- Nail diameter: Proximal 15.5mm, distal: 11.0mm
- Proximal Nail neck-shaft angles: 120°, 125°, 130°
- M-L bend for valgus curvature: 4 degrees
- Proximal anterversion: 10 degrees
- **End Caps:** Standard, +5mm and +10mm

#### **Long Nail Distal Locking Options**

**Dynamic locking** (Only one screw is needed):

• Locking in the distal part of the oblong hole creates a dynamic lokking mechanism (Figure 3).

#### Secondary dynamization

(Two screws are needed):

• One screw placed in the distal part of the oblong hole and the other in the round hole.

If dynamization is required after a period of time, the screw placed in the round hole has to be removed (Figure 4).

- **Radius of Curvature :** R2.0m of the shaft
- **Distal locking holes:** (round and oblong) for 5mm screws; up to 5mm of dynamization is possible
- Lag Screw diameter: 10.5mm
- Lag Screw lengths: 75-120mm in 5mm increments
- Patented Self-Tapping Lag Screw design for high load absorption and easy insertion
- **Asymmetrical depth profile** to allow the Lag Screw to slide in the lateral direction only for dynamic bone compression at the fracture site to enhance fracture healing
- A patented self retaining Set Screw to prevent rotation and simultaneously allowing sliding of the Lag Screw laterally
- Distal Locking Screw diameter: 5mm
- Fully threaded Distal Locking Screw lengths: ranging from 25-120mm in 2.5mm or 5mm increments
- Self-tapping screw tip with optimized short cutting flutes
- **Optimized diameter** under the head helps to prevent microfractures during insertion
- Low profile head

\* Each nail is supplied sterile, packaged together with a Set Screw in one box.

### Static locking (Two screws

are needed):

• One screw placed in the round hole and the other is placed in the proximal part of the oblong hole (Figure 5).







Gamma3 End Cap

Figure 3 Dynamic Locking

Figure 4 Secondary Dynamization

Figure 5 Static Locking

7

### Gamma3 Nail Indications

#### **Trochanteric Gamma3 Nail**

The Gamma3 Trochanteric Nail is indicated for:

- Intertrochanteric fractures
- Pertrochanteric fractures
- Nonunions and malunions

#### Contraindications for the Trochanteric Gamma3 Nail

- Medial neck fracture
- Subtrochanteric fracture



Cut-Away View

Lateral view of Lag Screw



Lag Screw Stabilization System

# Gamma3 Nail Indications

### Long Gamma3 Nail

The Long Gamma3 Nail is indicated for:

- Subtrochanteric fractures
- Pertrochanteric fractures associated with shaft fractures
- Pathological fractures (including prophylactic use) in both trochanteric and disphysal areas
- Nonunions and malunions

#### Contraindications for the Long Gamma3 Nail

• Medial neck fracture

### Long Gamma3 Nail Length Selection

The appropriate length of the Long Nail may be chosen either by pre-operative planning using an X-ray of the injured femur and a template, or intra-operatively as is the usual practice in intramedullary nailing.



#### **Implant Selection**

X-ray templates may be used for pre-operative planning. Alternatively the femoral neck angle may be determined through the use of a goniometer.



X-ray in A/P view, showing implant



(Cat. No 1320-0002)



Gamma3 Long Nail X-ray Template (Cat. No 1320-0005)

#### Note:

Ensure precise alignment of the affected hip joint when using these templates. Template magnification is 15%. All dimensions (nail angle and implant sizing) result from use of these templates and must be verified intraoperatively to ensure proper implant selection.

#### **Suggested Patient Positioning**

The patient is placed in a supine position on the fracture table and closed reduction of the fracture is recommended. Reduction should be achieved as anatomically as possible. If this is not achievable in a closed procedure, open reduction may be necessary.

Traction is applied to the fracture, keeping the leg straight. The unaffected leg is abducted to allow for use of the image intensifier (Figure 6).

While maintaining traction, the leg is internally rotated 10-15 degrees to complete fracture reduction; the patella should be horizontal or slightly inward.

Position the image intensifier so that A/P and M/L views of the trochanteric region of the affected femur may be easily obtained. This position is best achieved if the image intensifier is positioned so that the axis of rotation of the intensifier is centered on the femoral neck of the affected femur (Figure 7).

It is important to ensure that a view of both the distal and proximal ends of the nail can be obtained during the procedure without obstruction of the traction table.

The patient is then prepared and draped accordingly.



Figure 6



Figure 7

#### **Special Techniques for Fracture Reduction**

For specific situations, special techniques have been developed for fracture reduction, and these are explained below.

To counter this misalignment, the trunk is turned to the opposite side and held in position by a thoracic rest or by a large drape. This tightens the gluteus medius muscles and relaxes the psoas, externally rotating the proximal fragment into alignment and exposing the trochanter for easier introduction of the nail. The fractured limb is kept straight, with the knee in flexion (Figure 8), using the stirrup to avoid adduction. This position helps to align the distal portion. Reduction is confirmed in the A/P view.

Subtrochanteric fractures cannot always be reduced during positioning in the lateral view, because of the stresses applied to the proximal fragment by the psoas muscles. The fragment may be reduced during surgery by using the Universal Rod and the Reduction Spoon. (Figure 9).

Care must be taken when introducing the implant as the proximal fragment may rotate during insertion.







#### **Entry Point**

The correct entry point is located at the junction of the anterior third and posterior two-thirds of the tip of the greater trochanter and on the tip itself (Figure 10). A skin incision of appropriate length centered between the tip of the greater trochanter and the anterior superior iliac spine. Incision site is also determined by the amount of soft tissue around immediate area.

#### **Preparation of the Medullary Canal**

The Cannulated Curved Awl (Cat. No. 1806-0041, Figure 11) is used to open the medullary canal under image intensification. The optional Cannulated Cutter (Cat. No. 1320-0041) may be used in place of the Cannulated Curved Awl (Figure 12).



Figure 10



Figure 11



13

#### Trochanteric Gamma3 Nail Reaming the Medullary Canal

A 3mm Ball-Tipped Guide Wire (Cat. No. 1806-0085S) is passed through the Cannulated Curved Awl into the shaft of the femur as shown, using the Guide Wire Handle (Cat. No. 1806-0095 and 1806-0096, Figure 13).

Flexible reamers are used to ream the shaft of the femur in stages starting from 9mm diameter and increasing in 0.5mm increments (Figure 13a). The canal should be reamed at least 13mm (Figure 14).

In order to accommodate the proximal portion of the nail, the subtrochanteric region must be opened to 15.5mm (Figure 15). This can be done either by reaming with the Stryker Bixcut Reaming System (Figure 13a) or, alternatively, with the One Step Conical Reamer (Cat. No. 1320-0011). For soft tissue protection, the Conical Reamer Sleeve (Cat. No. 1320-0031) should be used during reaming.

Care must be taken with flexible reamers to ensure that the Guide Wire is not displaced laterally during reaming.



Figure 13a



Figure 14



#### Long Gamma3 Nail Reaming the Medullary Canal

A 3mm Ball-Tipped Guide Wire (Cat. No. 1806-0085S) is passed through the Cannulated Curved Awl into the shaft of the femur as shown, using the Guide Wire Handle (Cat. No. 1806-0095 and 1806-0096, Figure 16).

Flexible reamers are used to ream the shaft of the femur in stages starting from 9mm diameter and increasing in 0.5mm increments (Figure 16a). The canal should be reamed at least 2mm larger than the distal diameter of the nail (Figure 17).

When reaming is performed, the entire femoral canal should be over-reamed down through the isthmus, in order to reduce stress risers in the bone.

In order to accommodate the proximal portion of the nail, the subtrochanteric region must be opened to 15.5mm (Figure 18). This can be done either by reaming with the Stryker Bixcut Reaming System (Figure 16a) or, alternatively, with the One Step Conical Reamer (Cat. No. 1320-0011). For soft tissue protection, the Conical Reamer Sleeve (Cat. No. 1320-0031) should be used during reaming.

Care must be taken with flexible reamers to ensure that the Guide Wire is not displaced laterally during reaming.



Figure 16a



Figure 17



15

#### Using the Optional "One Step Conical Reamer"

The One Step Conical Reamer (Cat. No. 1320-0011) is an optional instrument that has been developed to provide surgeons with another option to prepare the proximal canal of the femur using only one drilling step.

When the Gamma3 Nail is used, reaming of the subtrochanteric and diaphyseal region of the femoral cavity may not be required, particularly in elderly patients with wide medullary canals.

After skin incision and positioning of the Guide Wire as described above, the Trocar (Cat. No. 1320-0021) or Multi Hole Trocar (Cat. No. 1320-0026) is inserted into the Conical Reamer Sleeve (Cat. No. 1320-0031) to protect the soft tissue during insertion. Push the Trocar (use center hole if Multi Hole Trocar is used) and Sleeve Assembly down over the 3mm Guide Wire to the tip of the trochanter (Figure 19 and 20).

### **Entry Point Optimization**

The Entry Point can also be made without using the Awl or the Cannulated Cutter. A 3.2mm K-Wire is placed through the tip of the trochanter.

If you find that the K-Wire is not positioned in the optimal position, it may be easily corrected using a second K-Wire in combination with the Multi Hole Trocar.

The Multi Hole Trocar (Figure 20a) has a special design for more precise insertion. In addition to the central hole, four other holes are located eccentrically at different distances from the center (Figure 20) to easily revise insertion of the guiding K-Wire into the proper position (Entry Point).



#### **Assembly of the Targeting Device**

#### 1. Targeting Sleeve and Knob Assembly

Align the arrow on the Targeting Sleeve (Cat. No. 1320-0118) with the "spot" on the Knob (Cat. No. 1320-0105) for assembly (Figure 21a). Slide and push the knob approximately 5mm onto the Targeting Sleeve and turn clockwise one-third of a turn to secure in place (Figure 21b).

#### 2. Targeting Arm and Targeting Sleeve Assembly

With the Knob in place on the Targeting Sleeve (Figure 21b) align the long arrow on the Targeting Sleeve with the long arrow on the Targeting Arm (Cat. No. 1320-1000) and assemble (Figure 21c). Rotate the Targeting Sleeve around the barrel of the Plus Targeting Arm until the appropriate nail angle (eg 125°, point-to-point) or distal locking position (static or dynamic) is reached. Once adjusted, the Targeting Sleeve is snapped into place by pushing towards the proximal part of the Targeting Device. Targeting Arm Cat. No. 1320-1000

Targeting Sleeve (green coded) Cat. No. 1320-0118





Gamma3<sup>™</sup> Plus



Knob Cat. No. 1320-0105

Figure 21b

Figure 21c



### 3. Assembly of the Targeting Device and the Gamma3 Nail

The appropriate size Gamma3 Nail is now assembled to the Targeting Device (Figure 22a). The Targeting Device is keyed and the nail should be aligned with the appropriate pegs for proper assembly.

Securely tighten the Nail Holding Bolt (Cat. No. 1320-0990) with the Ball Tip Screwdriver (Cat. No. 1320-0065), so that it does not loosen during nail insertion.

Before starting surgery the following functions of the Target Device should be confirmed:

- 1. There is a secure connection between the Nail and the Target Device.
- 2. The Lag Screw Guide Sleeve (Cat. No. 1320-0130) is correctly aligned with the selected nail angle (Figure 22b).
- 3. The Distal locking position of the Tissue Protection Sleeve (Cat. No. 1320-0215), is aligned for either static or dynamic locking (Figure 22c).



#### **Nail Insertion**

Insert the Gamma3 Nail by hand (Figure 23).

#### Optional Gamma3 Plus Final Seating Impactor

The carbon fiber guide should never be struck as it may break or become deformed. The impactor that is provided can be utilized to assist with final seating of the long nail (Figure 23a). Gentle tapping will produce small adjustments (in the nail position) that can help to optimize the ultimate position of the lag screw in the femoral head. Following any use of the impactor, the connecting bolt should be re-tightened.

The impactor should not be utilized to force the nail down the canal. If the nail cannot be seated manually or if there is no advancement each time the impactor is tapped, A/P and Lateral fluoroscopic X-rays should be reviewed to determine the cause of the impingement - there may be a mis-match between the nail geometry and the medullary canal. The starting position, the femoral bow and the canal diameter should all be examined to insure that the leading end of the nail is not impinging on the medial or anterior cortex and that the canal itself has been sufficiently reamed. Periodically, nail removal and further reaming of the diaphysis may be required.

In some cases (particularly young patients or those of short stature), the proximal metaphyseal flair may be undersized and serve to prevent nail advancement. If this situation is encountered, a 15.5 or 16mm flexible reamer may be used to further widen this area to (at least) the level of the lesser trochanter.

Beginning long nail insertion with the nail internally rotated 90° until the fracture has been passed may help to facilitate manual passage. The final nail depth position is monitored with the image intensifier. The Lag Screw position may be estimated with a ruler on the C-arm monitor to confirm proper Lag Screw position. The Lag Screw should be placed centrally or slightly inferior in the femoral head and frontal plane (Figure 24).

#### Note:

#### Remove the Guide Wire (Cat. No. 1806-0085S) following successful placement of the nail (Figure 25).

With the Gamma3 Nail in place, check the anteversion of the nail. If the K-Wire Clip (Cat. No. 1320-0110) is used to check anteversion, it is mounted onto the slots with the u wire (1320-0112) in the Targeting Arm by depressing the clip flanges (Figure 26).

#### Before proceeding ensure the Nail Holding Bolt is fully tightened to the nail.

The Lag Screw should be placed in the central position of the femoral head in the lateral view (Figure 27).



Figure 23



Figure 23a



Figure 24



Figure 25



Figure 26





#### **Optional Technique**

### Lag Screw Positioning using the One Shot Device

The One Shot Device is recommended for optimal Lag Screw placement:

The One Shot Device (Cat. No. 1320-3010) is recommended for establishing whether the Lag Screw is in the optimal position. This device enables correct positioning of the K-Wire for Lag Screw placement before performing the lateral skin incision and opening of the lateral cortex (Figure 28).



correct position



igure 28a A/P view









Figure 31 Lag Screw Guide Sleeve in good contact to the lateral cortex

#### **Lag Screw Insertion**

Insert the Lag Screw Guide Sleeve with the 4.2mm Drill Guide Sleeve (green, Cat No. 1320-0140) through the Targeting Sleeve to the level of the skin. Make the appropriate incision and advance the Guide Sleeve to the bone (Figure 29).

#### Note:

For accurate **Lag Screw measurement the Guide Sleeve must be in contact with the bone**. Once the Guide Sleeve is on the bone, it is locked in place to stabilize the Targeting Assembly, by turning the Knob clockwise (Figure 30 and 31).

#### **Lag Screw Insertion**

With the Lag Screw Guide Sleeve firmly in place the  $4.2 \times 300$  mm Center Tipped Drill (Cat No. 1320-3042S), is used to open the first cortex. (Figure 32).

The green coded 4.2mm Lag Screw Drill Guide Sleeve is then replaced by the K-Wire Sleeve (Cat. No. 1320-0150).

The sleeves look similar, but have different inner hole diameters. The K-Wire Sleeve has no colored ring (Figure 33).



The 3.2 x 450mm K-Wire (Cat No. 1210-6450S) is inserted through the K-Wire Sleeve and advanced to the bone (Figure 34), using the Guide Wire Handle (Cat No. 1806-0095 and 1806-0096). Check that the K-Wire is placed either central or in the lower half of the femoral head in the frontal plane and on the mid-line in the lateral plane (Figure 35).

Check the position with the image intensifier in both the A/P and M/L views to ensure proper K-Wire positioning.



Figure 32 Opening of the lateral cortex



Figure 35 K-Wire placement

![](_page_21_Picture_1.jpeg)

Figure 37 Lag Screw Length Measurement

#### **Lag Screw Insertion**

After proper positioning of the K-Wire, the required Lag Screw length is measured using the Lag Screw Ruler (Cat. No. 1320-0180).

#### Note:

Before measuring, ensure that the Lag Screw Guide Sleeve is still pressed firmly against the bone (Figure 36).

Measure the appropriate Lag Screw length directly off of the Guide Wire with the Lag Screw Ruler (Figure 37).

If the measurement falls between markings on the scale, e.g. 97mm, it should always be rounded up to the next higher value, e.g. 100mm.

![](_page_21_Picture_9.jpeg)

Figure 37a Suggested Pin Placement

Lock K-Wire window Figure 38

The Lag Screw length measurement is set with the adjustable stop on the Lag Screw Step Drill (Cat. No. 1320-0190). The appropriate length must be visible in the Step Drill "window". (Figure 38 & 39).

The K-Wire Sleeve is removed and the Lag Screw Step Drill, set to the appropriate length, is passed over the K-Wire, through the Lag Screw Guide Sleeve (Figure 40).

The channel for the Lag Screw may be prepared using the T-Handle (Cat. No. 702628) connected to the Lag Screw Step Drill.

Drilling should continue until the Step Drill Stop comes into contact with the Lag Screw Guide Sleeve (Figure 40a).

Drilling should be performed under image intensifier control to avoid penetration of the hip joint. The K-Wire may also be observed in the K-Wire window of the Step Drill (Figure 40b).

### Note:

It is important to observe the K-Wire tip during drilling on the C-Arm. Additionally, the K-Wire window provides another means for monitoring the K-Wire position.

![](_page_22_Picture_9.jpeg)

Figure 39

![](_page_22_Picture_11.jpeg)

![](_page_23_Picture_1.jpeg)

Check on image intensifier during drilling to monitor the depth of the drill as it approaches subchondral bone.

Approximately 8 to 10mm of K-Wire will remain outside of the Step Drill. This additional length of K-Wire is to provide a margin of safety so that the Step Drill does not penetrate the joint (Figure 41).

The appropriate length of the Lag Screw should be the same as that of the Step Drill. The screw is then assembled to the Lag Screwdriver (Cat. No. 1320-0200).

Figure 42 Lag Screw and Lag Screwdriver assembly

![](_page_23_Picture_7.jpeg)

![](_page_23_Picture_8.jpeg)

Figure 43a

Ensure that the Lag Screw is firmly attached to the Lag Screwdriver (Figure 42).

Compression across the fracture site may be applied by selecting a Lag Screw 5mm shorter than originally measured (see Compression/Apposition below).

The Lag Screw assembly is now passed over the K-Wire, through the Lag Screw Guide Sleeve and threaded into the femoral head. Check the end position of the Lag Screw on the image intensifier.

#### **Compression/Apposition**

If compression or apposition of the fracture gap is required, this can be achieved by gently turning the thumbwheel of the Lag Screwdriver clockwise against the Guide Sleeve (Figure 43). Before starting compression, make sure that the Lag Screw Guide Sleeve is unlocked to allow for free sliding. To unlock the Lag Screw Guide Sleeve, the Knob is turned counter-clockwise. In osteoporotic bone, care must be taken to prevent Lag Screw pullout in the femoral head. The Lag Screw chosen should be shorter depending on the expected amount of compression.

![](_page_24_Picture_1.jpeg)

#### **Lag Screw Fixation**

The handle of the Lag Screwdriver must either be parallel or perpendicular (90°) to the Targeting Arm so that the Set Screw will engage one of the grooves of the Lag Screw shaft (Figure 44).

If the T-Handle is not perpendicular or parallel to the Targeting Arm, turn it clockwise until it reaches this position.

### DO NOT TURN THE LAG SCREW COUNTER-CLOCKWISE.

If the K-Wire is inadvertently removed, the Lag Screw may still be inserted without it, provided that the Guide Sleeve is still in contact with the cortex.

#### Note:

The Lag Screw should be placed in the pre-drilled hole in order to provide maximum resistance to cut out. The Lag Screw should not be turned counter-clockwise after the final position is reached, otherwise the Lag Screw may lose bony contact at the tip.

The Set Screw is inserted through the Targeting Device and Nail Holding Bolt, and while in place, attached to the nail with the Set Screwdriver. Turn the Set Screw until contact is made with one of the grooves of the Lag Screw (Figure 45).

#### Note:

The handle of the Lag Screwdriver is either parallel or at right angles (90°) to the Targeting Arm for proper placement of the Set Screw.

To verify the correct position of the Set Screw, turn the Lag Screwdriver gently clockwise and counter-clockwise. If resistance is met, the Set Screw is in place. If the Lag Screw moves, re-align the handle and tighten the Set Screw again until it engages in one of the four grooves.

After slightly tightening the Set Screw it should then be backed out by one quarter (1/4) turn. This allows the Lag Screw to slide freely while maintaining rotational control of the Lag Screw (Figure 46).

#### Note:

Proper Set Screw placement is critical to the sliding characteristics of the Lag Screw. Do not unscrew the Set Screw more than ¼ turn.

If distal locking is not indicated, an End Cap may be inserted at this point to prevent bony ingrowth at this stage of the procedure. To insert the appropriate size End Cap the Lag Screwdriver is left in place as the Nail Holding Bolt is removed. Once the Nail Holding Bolt is removed an End Cap may be inserted through the Targeting Arm.

Alternatively the End Cap may be inserted free-hand after removal of the Targeting Device.

![](_page_25_Picture_7.jpeg)

Figure 46

#### **Distal Screw Locking**

Disconnect the Lag Screwdriver by loosening the end thumbwheel, remove the Lag Screwdriver, Lag Screw Guide Sleeve and the K-Wire. The nature of the fracture determines whether the distal Locking Screw is used.

It should be used:

- If the fracture is unstable
- If rotational stability is required
- When there is a wide disparity between the diameter of the nail and the femoral cavity

Gamma3 nails offer the possibility to be locked distally either dynamically or statically. The fracture pattern determines the method of distal locking.

#### Note:

### These following points must be considered in order to perform a proper distal locking procedure:

- Ensure that the Nail Holding Bolt is still fully tightened
- Avoid soft tissue pressure on the distal locking sleeve assembly; therefore, the skin incision should be made co-linear in the direction of the sleeve assembly
- Check that the distal locking sleeve assembly, with the trocar removed, is in contact with the lateral cortex of the Femur and is locked securely with the Knob. Confirm final locking screw placement with A/P and Lateral flouroscopic X-Ray views
- Neutralize the power tool weight during drilling procedure and do not apply force to the Targeting Arm
- Start the power tool before having bone contact with the drill
- Use sharp and center tipped drills only

The Carbon Fiber Targeting Device offers the options of guided distal locking in a dynamic or static position of the nail. The green coded Targeting Sleeve of the Target Arm has to be adjusted in the required position. In the following description, a dynamic locking will be described. Turn the Targeting Sleeve until you reach the dynamic position with the point on the Target Sleeve is in line where the arrow on the target arm. Push the sleeve up in the cranial direction. Now assemble the Distal Tissue Protector, Drill Guide Sleeve and Trocar and advance it through the hole of the Target Arm down to the skin. A small incision is started at the tip of the Trocar, and is extended down to the lateral cortex (Fig. 48). The Trocar will extend back of the sleeve by approx. 3mm when the Tissue Protection Sleeve has reached the lateral cortex (Fig. 49).

![](_page_26_Picture_19.jpeg)

![](_page_26_Picture_20.jpeg)

Before locking the sleeve, gently turn the Knob clockwise, making sure that the Tissue Protection Sleeve is in good contact to the bone (Fig. 49).

The Trocar is now removed and replaced by the calibrated green coded  $4.2 \text{mm} \times 300 \text{mm}$  drill. Drill through the first cortex and as the second cortex is reached read off the measurement on the drill scale. Add the thickness of the cortex, which is approximately 5mm, to this measurement to select the correct screw length (Fig. 49a).

Alternatively, the drill can be drilled through the second cortex and monitored by X-ray or image intensifier. The screw length can then be **read directly** from the scale on the drill (Fig. 49a). Proceed to drill the second cortex.

It is also possible to measure the correct screw length using the Screw Gauge after drilling through the second cortex. The Drill Guide Sleeve must be removed and the Screw Gauge may be advanced through the Tissue Protection Sleeve. Put the small hook behind the medial cortex and read the required locking screw length from the scale.

Insert the 5mm distal Locking Screw through the Distal Tissue Protector by using the 3.5mm Screwdriver until the mark on the Screwdriver shaft approaches the Protector; advance the screw head carefully until it is slightly in direct contact with the cortex (Fig. 49a).

#### Note:

When the mark on the Screwdriver shaft reaches the Tissue Protection Sleeve, this indicates that the screw head is near the cortex (Fig. 49a). Take care not to overscrew. The screw head should come just into contact with the cortex and resistance should be felt.

![](_page_27_Picture_8.jpeg)

![](_page_28_Picture_1.jpeg)

![](_page_28_Picture_2.jpeg)

Figure 51a

![](_page_28_Figure_3.jpeg)

→ Not in line with the Nail holes
→ This is the best position to drill

![](_page_28_Picture_5.jpeg)

### Distal Screw Locking of the Long Gamma3 Nail

For distal locking, the Long Gamma3 Nail offers the following three options (Figure 51).

- 1. **Dynamic locking** (Only one screw is needed):
  - Locking in the distal part of the oblong hole creates a dynamic lock-ing mechanism.

### 2. Secondary dynamization

(Two screws are needed):

- One screw is placed in the distal part of the oblong hole and the other in the round hole. If dynamization is required after a period of time, the screw placed in the round hole has to be removed.
- 3. **Static locking** (Two screws are needed):
  - One screw is placed in the round hole and the other is placed in the proximal portion of the oblong hole.

Disconnect the Lag Screwdriver by loosening the end thumbwheel, remove the Lag Screwdriver, Lag Screw Guide Sleeve and the K-Wire. The nature of the fracture determines whether the distal Locking Screw is used.

It should be used:

- If the fracture is unstable
- If rotational stability is required
- When there is a wide disparity between the diameter of the nail and the femoral cavity

### Visualizing the Distal Holes

The essential first step in distal targeting is to position the image intensifier so that the upper distal hole in the nail appears perfectly round. The oblong hole does not appear round. If the hole appears to be elliptical in either the vertical or horizontal planes, the image intensifier position must be adjusted appropriately as shown in (Figures 51a and 51b). It is advised to correct image in one plane at a time.

30

#### Distal Screw Locking of the Long Gamma3 Nail

#### **Free-hand Technique**

The free-hand technique is used to insert Locking Screws into the distal holes of the nail. Rotational alignment must be checked prior to locking the nail. Distal nail locking is described as follows, using the Static Locking mode according to Figure 51, Number 3.

The center-tipped 4.2 x 180mm Drill (Cat. No. 1806-4270S) is held at an oblique angle to the center of the locking hole (Figure 53). Upon X-ray verification, the Drill is placed perpendicular to the nail and drilled through the lateral and medial cortex. Confirm in both the A/P and M/L planes by X-ray that the Drill passes through the hole in the nail.

After drilling both cortices, the screw length may be read directly off of the Long Screw Scale (Cat. No. 1806-0365) at the green ring on the center tipped Drill.

It is also possible to measure the correct screw length using the Screw Gauge (Cat. No. 1806-0480) after drilling through the second cortex. Put the small hook behind the medial cortex and read the required locking screw length from the scale (Figure 53).

Insert the 5mm distal Locking Screw (Figure 54) using the 3.5mm Screwdriver and advance the screw carefully until it comes into contact with the cortex (Figure 55).

#### **Important Note:**

Take care not to over tighten. The screw head should come just into contact with the cortex and resistance felt.

![](_page_29_Figure_10.jpeg)

Figure 52

![](_page_29_Picture_12.jpeg)

Figure 53

![](_page_29_Figure_14.jpeg)

Figure 54

![](_page_29_Figure_16.jpeg)

Figure 55

#### **End Cap Insertion**

It is recommended to use an End Cap to prevent bony ingrowth.

Leave the Screwdriver for distal locking in place to stabilize the Targeting Device and remove the Nail Holding Bolt. Load the Standard Size End Cap (Cat. No. 3005-1100S) onto the Screwdriver (Cat. No. 1320-0065) and pass the assembly through the top of the Targeting Device down into the nail (Figure 56).

Alternatively, the End Cap may be inserted free hand after removal of the Targeting Device.

End Cap (size Standard) Cat. No. 3005-1100S

![](_page_30_Picture_6.jpeg)

Once the End Cap is in place, remove the Distal Screwdriver and Targeting Device (Figure 57). The wounds are closed in the usual manner.

![](_page_30_Picture_8.jpeg)

Figure 56 End Cap assembly

![](_page_30_Picture_10.jpeg)

Figure 57 Final Nail assembly

#### **Nail Extension End Caps**

If the proximal end of the nail is counter sunk into the bone and cortical bone support is required, End Caps in sizes +5mm and +10mm are available.

![](_page_31_Picture_3.jpeg)

Extended End Caps may only be inserted free hand after the Targeting Device has been removed.

![](_page_31_Picture_5.jpeg)

End Cap (size +5mm) (Cat. No. 3005-1105S)

![](_page_31_Picture_7.jpeg)

End Cap (size Standard) (Cat. No. 3005-1100S)

Set Screw (Cat. No. 3003-0822S)

![](_page_31_Picture_10.jpeg)

End Cap (size +10mm) (Cat. No. 3005-1110S)

![](_page_31_Picture_12.jpeg)

![](_page_31_Picture_13.jpeg)

End Cap (size +5mm) (Cat. No. 3005-1105S)

End Cap (size Standard)

(Cat. No. 3005-1100S)

Set Screw (Cat. No. 3003-0822S)

![](_page_31_Figure_17.jpeg)

![](_page_31_Picture_18.jpeg)

![](_page_32_Picture_1.jpeg)

Figure 58

![](_page_32_Picture_3.jpeg)

Figure 59

![](_page_32_Picture_5.jpeg)

Figure 60

### Extraction of the Gamma3 Implants

When implant extraction is indicated, please proceed as follows:

#### Step I (Figure 58)

Remove the distal screw using the 3.5mm Screwdriver after making an incision through the old scar.

#### Step II (Figure 59)

Make a small incision through the old scar below the greater trochanter to expose the outer end of the Lag Screw. Remove any bony ingrowth, which may be obstructing the outer end or internal thread of the Lag Screw.

The K-Wire is then introduced via the Lag Screw into the head of the femur. The Lag Screwdriver is passed over the K-Wire, using the Lag Screw Guide Sleeve as a Tissue Protector, and engaged with the distal end of the Lag Screw.

Check that ingrowth does not obstruct secure engagement of the Lag Screwdriver, otherwise the Lag Screw or Screwdriver may be damaged and extraction may be much more difficult. Tighten the Lag Screwdriver to the Lag Screw by turning the thumbwheel clockwise.

#### **Step III (Figure 60)**

An incision is made over the proximal end of the nail. The proximal End Cap if present is removed using the Ball Tip Screwdriver or Strike Plate (Cat. No. 1320-0070), and the Set Screwdriver is engaged into the Set Screw. The screw is rotated counterclockwise until it is removed.

#### Step IV (Figure 61)

The Conical Extraction Rod is threaded and tightened into the proximal end of the nail. The Lag Screw is extracted by counter-clockwise rotation and pulling of the Lag Screwdriver. The K-Wire must then be removed.

#### Step V (Figure 62)

An appropriate Sliding Hammer assembly is attached to the Extraction Rod and the nail is extracted.

#### Note:

It may be useful to turn the Lag Screw Screwdriver slightly clockwise to loosen any bony ingrowth into the screw threads before turning it counter-clockwise.

![](_page_33_Picture_7.jpeg)

Figure 61

![](_page_33_Picture_9.jpeg)

![](_page_33_Picture_10.jpeg)

Figure 62

# Ordering Information **Titanium Implants**

目

### Gamma3 Trochanteric Nail Kit, Ti (includes Set Screw)

| Cat. No.   | Diameter<br>mm | Length<br>mm |  |
|------------|----------------|--------------|--|
| 3120-1180S | 15.5/11        | 180mm × 120° |  |
| 3125-1180S | 15.5/11        | 180mm × 125° |  |
| 3130-1180S | 15.5/11        | 180mm × 130° |  |
|            |                |              |  |

![](_page_34_Picture_3.jpeg)

| Cat. No.   | Diameter<br>mm | Length<br>mm |
|------------|----------------|--------------|
| 0000 00750 | 10 5           | 75-20-20     |
| 3060-00755 | 10.5           | 73000        |
| 3060-0080S | 10.5           | 80mm         |
| 3060-0085S | 10.5           | 85mm         |
| 3060-0090S | 10.5           | 90mm         |
| 3060-0095S | 10.5           | 95mm         |
| 3060-0100S | 10.5           | 100mm        |
| 3060-0105S | 10.5           | 105mm        |
| 3060-0110S | 10.5           | 110mm        |
| 3060-0115S | 10.5           | 115mm        |
| 3060-0120S | 10.5           | 120mm        |

### Set Screw, Ti (available separately)

| Cat. No.   | Diameter<br>mm | Length<br>mm | 間等 |
|------------|----------------|--------------|----|
| 3003-0822S | 8              | 17.5mm       |    |

### 5mm Fully Threaded Locking Screw, Ti

| Cat. No.   | Diameter<br>mm | Length<br>mm |
|------------|----------------|--------------|
| 1896-5025S | 5.0            | 25           |
| 1896-5027S | 5.0            | 27.5         |
| 1896-5030S | 5.0            | 30           |
| 1896-5032S | 5.0            | 32.5         |
| 1896-5035S | 5.0            | 35           |
| 1896-5037S | 5.0            | 37.5         |
| 1896-5040S | 5.0            | 40           |
| 1896-5042S | 5.0            | 42.5         |
| 1896-5045S | 5.0            | 45           |
| 1896-5050S | 5.0            | 50           |
| 1896-5055S | 5.0            | 55           |
| 1896-5060S | 5.0            | 60           |
| 1896-5065S | 5.0            | 65           |
| 1896-5070S | 5.0            | 70           |
| 1896-5075S | 5.0            | 75           |
| 1896-5080S | 5.0            | 80           |
| 1896-5085S | 5.0            | 85           |
| 1896-5090S | 5.0            | 90           |
| 1896-5095S | 5.0            | 95           |
| 1896-5100S | 5.0            | 100          |
| 1896-5105S | 5.0            | 105          |
| 1896-5110S | 5.0            | 110          |
| 1896-5115S | 5.0            | 115          |
| 1896-5120S | 5.0            | 120          |

### End Caps, Ti

|   | Cat. No.                               | Diameter<br>mm       | Length<br>mm          |
|---|--|----------------------|-----------------------|
| 1 | 3005-1100S<br>3005-1105S<br>3005-1110S | 11.0<br>15.5<br>15.5 | Standard<br>+5<br>+10 |
|   |  |                      |                       |

# Ordering Information **Titanium Implants**

#### Long Gamma3 Nail Kit R2.0, Ti (includes Set Screw) Right

| Cat. No.    | Diameter<br>mm | Length<br>mm           |  |
|-------------|----------------|------------------------|--|
| 3220-0280S  | 11             | 280 x 120°, Right      |  |
| 3220-0300S  | 11             | 300 x 120°, Right      |  |
| 3220-0320S  | 11             | 320 x 120°, Right      |  |
| 3220-0340S  | 11             | 340 x 120°, Right      |  |
| 3220-0360S  | 11             | 360 x 120°, Right      |  |
| 3220-0380S  | 11             | 380 x 120°, Right      |  |
| 3220-0400S  | 11             | 400 x 120°, Right      |  |
| 3220-0420S  | 11             | 420 x 120°, Right      |  |
| 3220-0440S  | 11             | 440 x 120°, Right      |  |
| 3220-0460S  | 11             | 460 x 120°, Right      |  |
| 3220-0480S* | 11             | 480 x 120°, Right      |  |
| 3225-02805  | 11             | 280 x 125° Bight       |  |
| 3225-0300S  | 11             | 300 x 125°, Right      |  |
| 3225-0320S  | 11             | 320 x 125°. Right      |  |
| 3225-0340S  | 11             | 340 x 125°. Right      |  |
| 3225-0360S  | 11             | 360 x 125°, Right      |  |
| 3225-0380S  | 11             | 380 x 125°, Right      |  |
| 3225-0400S  | 11             | 400 x 125°, Right      |  |
| 3225-0420S  | 11             | 420 x 125°, Right      |  |
| 3225-0440S  | 11             | -<br>440 x 125°, Right |  |
| 3225-0460S  | 11             | 460 x 125°, Right      |  |
| 3225-0480S* | 11             | 480 x 125°, Right      |  |
| 3230-02805  | 11             | 280 v 130° Right       |  |
| 3230-0300S  | 11             | 300 x 130° Right       |  |
| 3230-03205  | 11             | 320 x 130° Right       |  |
| 3230-0340S  | 11             | 340 x 130° Right       |  |
| 3230-0360S  | 11             | 360 x 130°. Right      |  |
| 3230-0380S  | 11             | 380 x 130°. Right      |  |
| 3230-0400S  | 11             | 400 x 130°. Riaht      |  |
| 3230-0420S  | 11             | 420 x 130°, Right      |  |
| 3230-0440S  | 11             | 440 x 130°, Riaht      |  |
| 3230-0460S  | 11             | 460 x 130°, Right      |  |
| 3230-0480S* | 11             | 480 x 130°, Right      |  |

#### Long Gamma3 Nail Kit R2.0, Ti (includes Set Screw) Left

| Cat. No.    | Diameter<br>mm | Length<br>mm     |
|-------------|----------------|------------------|
| 3320-0280S  | 11             | 280 x 120°, Left |
| 3320-0300S  | 11             | 300 x 120°, Left |
| 3320-0320S  | 11             | 320 x 120°, Left |
| 3320-0340S  | 11             | 340 x 120°, Left |
| 3320-0360S  | 11             | 360 x 120°, Left |
| 3320-0380S  | 11             | 380 x 120°, Left |
| 3320-0400S  | 11             | 400 x 120°, Left |
| 3320-0420S  | 11             | 420 x 120°, Left |
| 3320-0440S  | 11             | 440 x 120°, Left |
| 3320-0460S  | 11             | 460 x 120°, Left |
| 3320-0480S* | 11             | 480 x 120°, Left |
| 3325-0280S  | 11             | 280 x 125°, Left |
| 3325-0300S  | 11             | 300 x 125°, Left |
| 3325-0320S  | 11             | 320 x 125°, Left |
| 3325-0340S  | 11             | 340 x 125°, Left |
| 3325-0360S  | 11             | 360 x 125°, Left |
| 3325-0380S  | 11             | 380 x 125°, Left |
| 3325-0400S  | 11             | 400 x 125°, Left |
| 3325-0420S  | 11             | 420 x 125°, Left |
| 3325-0440S  | 11             | 440 x 125°, Left |
| 3325-0460S  | 11             | 460 x 125°, Left |
| 3325-0480S* | 11             | 480 x 125°, Left |
| 3330-0280S  | 11             | 280 x 130°, Left |
| 3330-0300S  | 11             | 300 x 130°, Left |
| 3330-0320S  | 11             | 320 x 130°, Left |
| 3330-0340S  | 11             | 340 x 130°, Left |
| 3330-0360S  | 11             | 360 x 130°, Left |
| 3330-0380S  | 11             | 380 x 130°, Left |
| 3330-0400S  | 11             | 400 x 130°, Left |
| 3330-0420S  | 11             | 420 x 130°, Left |
| 3330-0440S  | 11             | 440 x 130°, Left |
| 3330-0460S  | 11             | 460 x 130°, Left |
| 3330-0480S* | 11             | 480 x 130°, Left |

\*Special order only

### **Gamma3 Plus Instruments**

![](_page_36_Figure_2.jpeg)

### **Gamma3 Plus Instruments**

|           | Cat. No    | Description  |
|-----------|------------|--|
|           | 1806-0041  | Awl, Curved, Cannulated                                  |
|           | 1806-0085S | Guide Wire, Ball Tipped 3 x 1000mm, Sterile              |
|           | 1806-0095  | Guide Wire Handle  |
|           | 1806-0096  | Guide Wire Handle Chuck                                  |
|           | 1806-0185  | Tissue Protection Sleeve, Long (Distal)                  |
|           | 1320-0215  | Drill Sleeve, Long (Distal)                              |
|           | 1806-0232  | Screwdriver, Long  |
| <u></u> ا | 1320-0315  | Trocar, Long (Distal)                                    |
|           | 1806-0325  | Screw Gauge, Long  |
|           | 1806-0365  | Screw Scale, Long (for Long Nail)                        |
|           | 1806-0480  | Screw Gauge (for Long Nail)                              |
|           | 1806-4270S | Drill, 4.2 x 180mm, AO Small, Sterile<br>(for Long Nail) |

| 1320-9210 | Instrument Set, Plus, Complete with Instruments |
|-----------|---|
| 1320-9200 | Instrument Tray, Basic, Empty                   |
| 1320-0002 | X-Ray Template Trochanteric Gamma3              |
| 1320-0005 | X-Ray Template Long Gamma3, R 2.0               |

### **Gamma3 Optional Instruments**

|          | Cat. No                | Description   |
|----------|------------------------|---|
| <u> </u> | 1213-90918             | Guide Pin for Cannulated Cutter, 4.0 x 400mm, Sterile                     |
|          | 1320-0011              | One Step Conical Reamer   |
|          | 1320-0021<br>1320-0022 | Conical Reamer, Trocar, Short (Shown)<br>Conical Reamer, Trocar, Long     |
|          | 1320-0026<br>1320-0027 | Multi-hole Trocar, Short (Shown)<br>Multi-hole Trocar, Long               |
|          | 1320-0031<br>1320-0032 | Conical Reamer Sleeve, Short (Shown)<br>Conical Reamer Sleeve, Long       |
|          | 1320-0041              | Cannulated Cutter   |
| 0]       | 1320-0042              | Sleeve for Cannulated Cutter  |
| £        | 1320-0070              | Screwdriver Strike Plate  |
|          | 1320-0080              | Universal Socket Wrench (for Nail Holding Bolt)                           |
|          | 1320-0131              | Lag Screw Guide Sleeve, Navigated   |
| 2000     | 1320-0135              | Adaptor for One Shot Device, Gamma  |
|          | 1320-0160              | Fragment Control Clip<br>not compatible with Gamma3 Plus Targeting System |
| 4        | 1320-0170              | Fragment Control Sleeve   |
|          | 1320-3010              | One Shot Device, Gamma3   |
|          | 1320-30308             | Drill, 3.0 x 300mm, AO Small, Sterile                                     |
| €¢₽      | 1407-4006              | Nail Extraction Adapter   |

### **Gamma3 Optional Instruments**

|            | Cat. No    | Description  |
|------------|------------|--|
|            | 1806-0020  | Guide Wire Ruler (for Long Nail)                         |
|            | 1806-0110  | Universal Rod  |
|            | 1806-0125  | Reduction Spoon  |
| 5 <u> </u> | 1806-0130  | Wrench 8mm/10mm  |
|            | 1806-0170  | Slotted Hammer   |
|            | 1806-0450  | Tissue Protection Sleeve                                 |
|            | 1806-0460  | Drill Sleeve 4.2mm                                       |
|            | 1806-4290S | Drill, 4.2 x 230mm, AO Small, Sterile<br>(for Long Nail) |
|            | 702634     | Large A/O Hall Adaptor                                   |
|            |            |  |

1320-9005

Optional Instrument Tray, Empty

![](_page_40_Picture_0.jpeg)

| <br> |
|------|
|      |
|      |
| <br> |
|      |
|      |
| <br> |
|      |
|      |
| <br> |
|      |
|      |
|      |
|      |
|      |
|      |
|      |
|      |
|      |
|      |
|      |
|      |
|      |
|      |
|      |
|      |
|      |
|      |
|      |
|      |
|      |
|      |
|      |
|      |
|      |
|      |
|      |
|      |
|      |
|      |
|      |
|      |

| Gamma3 Nail<br>Notes |
|----------------------|
|                      |
|                      |
|                      |
|                      |
|                      |
|                      |
|                      |
|                      |
|                      |
|                      |
|                      |
|                      |
|                      |
|                      |
|                      |
|                      |
|                      |
|                      |
|                      |
|                      |
|                      |
|                      |
|                      |
|                      |

![](_page_42_Picture_0.jpeg)

Joint Replacements

Trauma, Extremities & Deformities

Craniomaxillofacial

Spine

**Biologics** 

**Surgical Products** 

Neuro & ENT

**Interventional Pain** 

Navigation

Endoscopy

Communications

Imaging

**Patient Handling Equipment** 

**EMS Equipment** 

325 Corporate Drive Mahwah, NJ 07430 **t: 201 831 5000** 

www.stryker.com

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Gamma3 and Stryker. All other trademarks are trademarks of their respective owners or holders.

Literature Number: LG3P-OT Rev. 1 MS/GS 02/09

Copyright © 2009 Stryker Printed in USA