



ALLOCLASSIC[®] **HIP SYSTEM**

Surgical Technique

Classic proven design
with superior clinical
results.



zimmer
Confidence in your hands™

ALLOCLASSIC HIP SYSTEM

Developed in conjunction with:

Frank Gottschalk, MD

Professor, Orthopedic Surgery
The University of Texas
Southwestern Medical Center at Dallas
Dallas, TX

D. Kevin Lester, MD

Assistant Professor, Orthopedic Surgery
University of California at San Francisco
Fresno, CA

CONTENTS

| | |
|--|----|
| INTRODUCTION | 3 |
| PREOPERATIVE TEMPLATING/PLANNING | 4 |
| POSITIONING FOR X-RAYS | 4 |
| TEMPLATING THE FEMUR | 4 |
| FEMORAL PREPARATION | 6 |
| OSTEOTOMY OF THE FEMORAL NECK | 6 |
| PREPARATION OF THE FEMORAL CANAL | 7 |
| TRIAL REDUCTION | 11 |
| TRIAL OPTION | 12 |
| STEM IMPLANTATION | 13 |
| IMPLANT SPECIFICATIONS AND DIMENSIONS .. | 14 |
| ORDERING INFORMATION | 16 |



ALLOCLASSIC® HIP SYSTEM

The Alloclassic stem was introduced to the global market in 1986, after seven years of extensive research and development. Since its introduction, the Alloclassic has established high standards for total hip arthroplasty with an unchanged design, excellent clinical results, and a technically simple application.

Proven Classic

The long-term clinical success of the Alloclassic can be attributed to a unique combination of design features. A tapered stem geometry and a grit-blasted surface provide proven initial and secondary fixation. In addition, the Alloclassic's rectangular cross-section has unparalleled rotational stability while the broaching technique allows for minimal disturbance to the endosteal blood supply and maximum bone conservation.

Excellent Long-Term Clinical Results

The Alloclassic Hip is regarded by surgeons world-wide as the standard for total hip arthroplasty with over 16 years of excellent clinical results.

Broad Product Indications

The versatile Alloclassic can be used in trauma cases, primary hip arthroplasty and revision hip arthroplasty in all bone types, from low- to high-demand patients.

Compatibility

The Alloclassic's 12/14 taper is compatible with multiple Zimmer acetabular options, including unipolar and bipolar heads, ceramic heads, Metasul® metal-on-metal and Durasul® highly crosslinked polyethylene in head diameters up to 44 mm.

Biologic Fixation

The Alloclassic's grit-blasted titanium niobium alloy with a 4-8 micron surface roughness provides a proven surface for reliable biologic fixation for long-term fixation of the implant.

MIS Compatibility

The Alloclassic's tapered design and non-reaming technique allow for minimal disturbance to the endosteal bloody supply and maximum bone conservation, and make it the perfect choice for Minimally Invasive Surgery. The Alloclassic hip has been successfully used in MIS hip arthroplasty for over 16 years.

Extensive Sizing

Fourteen anatomical sizes allow for an optimal patient fit.

Preoperative Planning

The primary objectives of preoperative planning are to:

1. Determine preoperative leg length discrepancy.
2. Assess acetabular component size and placement.
3. Determine femoral component size, position and fit.
4. Assess the necessary femoral offset.

In addition, preoperative planning will assist in the identification of bone abnormalities and potential problems before surgery that might require special instrumentation.


Positioning for X-rays

For the A/P X-ray of the pelvis, the femurs should be internally rotated 15° to show an accurate view of the femoral neck length, metaphysis and diaphysis. A direct lateral X-ray may also be beneficial in determining implant sizing. The Alloclassic Hip System templates incorporate 18 percent magnification.

Templating the Femur

To estimate the leg length discrepancy on the X-ray a line should be drawn through the bottom of the ischium (Fig. 1). The distance should then be measured from the lesser trochanter to the drawn reference line. The measured difference between each measured side is the radiographic leg length discrepancy. As an alternate reference point, the tip of the greater trochanter to the drawn reference line may be measured.

Fig. 1



Templating the Femur

Select the femoral template size that will best fit the proximal and distal femur and equalize the leg length. The tapered geometry of the Alloclassic should fill the canal from the medial to the lateral cortical wall (Fig. 2). The femoral template should be in line with the long axis of the femur, in the neutral position. Draw the neck resection line at the point where the selected stem provides the desired amount of leg length. The proximal tip of the prosthesis and the tip of the greater trochanter are suitable reference points for determining the height of the final implant.

Fig. 2



Preoperative A/P View



Postoperative A/P View



Postoperative Lateral View

Femoral Preparation

Osteotomy of the Femoral Neck

The osteotomy of the femoral neck is 1-2 cm above the lesser trochanter at an angle of 40-45° to the neutral axis of the femur or parallel to the intertrochanteric line. This may vary due to differences in the proximal femoral anatomy and should be based on preoperative planning (Fig. 3-4).

The femoral head is then removed (Fig. 5).

The leg is then externally rotated to neutral in full extension to provide the best exposure for preparation of the acetabulum.

Fig. 3




Fig. 4



Fig. 5



Instruments Used

Box Chisel [01.00129.650]



Preparation of the Femoral Canal

After the acetabular component has been placed, the leg is then repositioned for optimal access to the femoral canal. With the knee bent at 90°, the leg is internally rotated. The cut surface of the femoral neck is now easily accessible for femoral preparation, for a posterior approach.

With the antero-lateral approach, the leg is adducted and externally rotated to allow for access to the cut surface of the femoral neck.

A rectangular box chisel is used to cut a slot in the proximal neck and trochanteric region.

The chisel should be as lateral as possible against the transition to the greater trochanter to create a slot to accommodate the prosthesis in a neutral position (Fig. 6-7).

Fig. 6



Fig. 7



Preparation of the Femoral Canal

The slot may be extended laterally by notching the cortex of the trochanter using a rongeur or small rasp (Fig. 8 and 9).

Next, the canal finder should be used to further open the medullary canal. Caution should be used to ensure that a neutral opening to the canal is established and varus positioning is avoided (Fig. 10).


Fig. 8



Fig. 9



Fig. 10



Instruments Used

Adaptor with Lever [6835]



Preparation of the Femoral Canal

Prepare the femoral canal by first using rasp size 1, and progress to the next larger size rasp. Only in cases of very small femora should one start with the special sizes 01 and 0. To ensure proper rasp height and subsequent implant height, calibration marks are located on the connector piece, which reference the height of the final implant (Fig. 11).


The proximal tip of the Alloclassic will rest just below the tip of the greater trochanter when the implant is finally seated.

The shaft of the slaphammer is aligned along the femoral axis, while the impactor handle may be used to reference the anteversion of the rasp and final hip implant (Fig. 12).

Fig. 11



Fig. 12



Preparation of the Femoral Canal

The insertion of the first rasp will determine the anteversion of the subsequent rasps and the final implant, so precise orientation is necessary for exact stem positioning (Fig. 13).

Once again, the lateral trochanter area should be periodically rasped to avoid varus implantation. Each rasp should be impacted to the level of the osteotomy and the angled proximal surface of the rasp.

The first contact with the cortex usually occurs at the distal end of the rasp. Progress will become more difficult with each progressive impaction and the pitch of the impactor will change from low to high once cortical bone is engaged. The next larger size rasp is usually the final size as preoperatively templated. The final rasp should be seated to the resection line and there should be no visible movement of the broach when the slaphammer is rotated.

The final implant size will directly correspond to the final rasp size.

Fig. 13



Instruments Used

Tapered Bolt [01.00209.014]




Trial Reduction

Trial reduction of the Alloclassic may be accomplished with either the rasp or the trial prosthesis. The two methods are slightly different, so they are each described below. The recommended method for trialing the Alloclassic is with the trial prosthesis because this best replicates the final implant range of motion and stem positioning.

Trialing directly off of the rasp is first accomplished by removing the rasp handle and leaving the rasp in the femoral canal. Either a trial rasp ball head or a trial rasp taper is then inserted into the rasp hole (Fig. 14). If the trial rasp taper is selected, any of the standard 12/14 trial ball heads may then be used in conjunction with a trial rasp taper. A special trial rasp taper is available for the smallest rasps, sizes 01 and 0 (Part #6842).

Once the appropriate trial head is selected, the hip is reduced. Leg length and offset are checked. This procedure is repeated as necessary using different length trial heads until optimal offset and leg length are established. A trial reduction should not allow significant push-pull or “shuck” of the joint in full extension. Range of motion is checked to avoid bony impingement and instability.

Fig. 14



Trial Option

Trialing with the trial prosthesis is first accomplished by removing the final rasp. A trial prosthesis matching the size of last rasp is inserted and driven in with the impactor instrument on the proximal tip of the trial prosthesis (Fig. 15). Once the trial is firmly seated, the trial ball head size as templated is connected to the trial taper (Fig. 16). The hip is then reduced. Next, leg length and offset are checked. This procedure is repeated as necessary using different length trial heads until optimal offset and leg length are established. A trial reduction should not allow significant push-pull or “shuck” of the joint in full extension. Range of motion is checked to avoid bony impingement and instability. The trial prosthesis is removed with the slaphammer and extractor hook, which is inserted into the proximal tip of the trial prosthesis (Fig. 17).

Fig. 15



Fig. 16



Fig. 17



Stem Implantation

The stem is inserted and impacted using the impactor in contact with the proximal tip of the prosthesis (Fig. 18).

The taper protector is removed from the taper and a trial head may be applied at this time for a final trial head reduction. Once the final range of motion and shuck tests are complete, the taper is carefully cleaned. The selected cobalt chrome or ceramic 12/14 taper femoral head is mounted with a light rotational movement, and rotated further with axial force until it is firmly seated. The ball head is seated with several taps with the ball head impactor instrument (Fig. 19).


Fig. 18



Fig. 19



Implant Specifications and Dimensions




| Size | Implant Length (mm) | Stem Length (mm) | Neck Length (mm) | Offset (mm) | Neck Angle | Proximal Width (mm) | | Distal Width (mm) | | Distal Section Distance (mm) |
|------|---------------------|------------------|------------------|-------------|------------|---------------------|----------|-------------------|----------|------------------------------|
| | | | | | | A-P View | M-L View | A-P View | M-L View | |
| 01 | 130 | 110 | 24 | 33 | 131° | 31 | 9 | 7 | 6 | 119 |
| 0 | 134 | 113 | 25 | 34 | 131° | 32 | 10 | 8 | 6 | 122 |
| 1 | 138 | 117 | 26 | 35 | 131° | 32 | 10 | 9 | 7 | 125 |
| 2 | 143 | 121 | 26 | 36 | 131° | 33 | 11 | 10 | 7 | 127 |
| 3 | 147 | 125 | 27 | 37 | 131° | 34 | 11 | 10 | 7 | 130 |
| 4 | 151 | 129 | 28 | 38 | 131° | 36 | 11 | 11 | 8 | 133 |
| 5 | 156 | 133 | 29 | 40 | 131° | 37 | 12 | 12 | 8 | 136 |
| 6 | 161 | 137 | 30 | 41 | 131° | 38 | 12 | 13 | 8 | 139 |
| 7 | 166 | 142 | 31 | 42 | 131° | 40 | 12 | 14 | 9 | 142 |
| 8 | 172 | 147 | 32 | 44 | 131° | 41 | 13 | 16 | 9 | 145 |
| 9 | 178 | 152 | 34 | 46 | 131° | 42 | 14 | 17 | 10 | 148 |
| 10 | 184 | 157 | 35 | 47 | 131° | 44 | 14 | 18 | 10 | 151 |
| 11 | 190 | 163 | 36 | 49 | 131° | 46 | 15 | 20 | 11 | 155 |
| 12 | 197 | 168 | 37 | 50 | 131° | 48 | 16 | 22 | 12 | 159 |

Implant Characteristics

The Alloclassic hip is the most widely used primary cementless hip in the world. Over 25,000 are implanted each year with over 350,000 Alloclassics implanted since the 1986 market introduction.

- Straight, collarless, tapered wedge design
- Rectangular cross section
- Grit blasted PROTASUL™-100 (Ti-6Al-7Nb)
- Surface Roughness $R_a = 251\mu\text{in.}$ or $(6.38\mu\text{m.})$
- Standard 12/14 taper

Implant Specifications and Dimensions



| Size | Implant Length (mm) | Stem Length (mm) | Neck Length (mm) | Offset (mm) | Neck Angle | Proximal Width (mm) | | Distal Width (mm) | | Distal Section Distance (mm) |
|------|---------------------|------------------|------------------|-------------|------------|---------------------|----------|-------------------|----------|------------------------------|
| | | | | | | A-P View | M-L View | A-P View | M-L View | |
| 3 | 178 | 156 | 31 | 40 | 131° | 34 | 11 | 10 | 7 | 130 |
| 4 | 183 | 160 | 32 | 41 | 131° | 36 | 11 | 11 | 7 | 133 |
| 5 | 188 | 165 | 33 | 43 | 131° | 37 | 12 | 12 | 8 | 136 |
| 6 | 193 | 169 | 34 | 44 | 131° | 38 | 12 | 13 | 8 | 139 |
| 7 | 199 | 174 | 35 | 45 | 131° | 40 | 12 | 14 | 9 | 142 |
| 8 | 204 | 179 | 36 | 47 | 131° | 41 | 13 | 16 | 9 | 145 |
| 9 | 210 | 184 | 38 | 48 | 131° | 42 | 14 | 17 | 10 | 148 |
| 10 | 215 | 189 | 39 | 50 | 131° | 44 | 14 | 18 | 10 | 151 |
| 11 | 221 | 194 | 40 | 52 | 131° | 46 | 15 | 20 | 11 | 155 |

Implant Characteristics

The Alloclassic SLL revision hip has 20% more distal length and is 4 mm longer in the neck than the Alloclassic SL primary hip.

- Straight, collarless, tapered wedge design
- Rectangular cross section
- Grit blasted PROTASUL™-100 (Ti-6Al-7Nb)
- Surface Roughness $R_a = 251\mu\text{in.}$ or $(6.38\mu\text{m.})$
- Standard 12/14 taper

Ordering Information

ALLOCLASSIC HIP IMPLANTS

| Catalog No. | Description |
|-------------|------------------------------|
| 2839 | Alloclassic SL Hip - Size 01 |
| 2840 | Alloclassic SL Hip - Size 0 |
| 2841 | Alloclassic SL - Size 1 |
| 2842 | Alloclassic SL - Size 2 |
| 2843 | Alloclassic SL - Size 3 |
| 2844 | Alloclassic SL - Size 4 |
| 2845 | Alloclassic SL - Size 5 |
| 2846 | Alloclassic SL - Size 6 |
| 2847 | Alloclassic SL - Size 7 |
| 2848 | Alloclassic SL - Size 8 |
| 2849 | Alloclassic SL - Size 9 |
| 2850 | Alloclassic SL - Size 10 |
| 2851 | Alloclassic SL - Size 11 |
| 2852 | Alloclassic SL - Size 12 |

ALLOCLASSIC SL HIP INSTRUMENT CASE 1

| Catalog No. | Description |
|--------------|--|
| 01.00129.600 | Instrument Case |
| 01.00129.601 | Insert for Instrument Case |
| 01.00029.031 | Instrument for Case Cover |
| 8095 | Striking Platform/Handle, Alloclassic |
| 6842 | Tapered Bolt (sizes 01 & 0), Alloclassic |
| 01.00209.014 | Tapered Bolt (sizes 1-12), Alloclassic |
| 5112 | Extractor Hook |
| 5671 | Gliding Hammer |
| 5672 | Impactor |
| 01.00129.650 | Chisel |
| 6835 | Adapter with Lever |
| 6815 | Detachable Rasp - Size 01 |
| 6820 | Detachable Rasp - Size 0 |
| 6821 | Detachable Rasp - Size 1 |
| 6822 | Detachable Rasp - Size 2 |
| 6823 | Detachable Rasp - Size 3 |
| 6824 | Detachable Rasp - Size 4 |
| 6825 | Detachable Rasp - Size 5 |
| 6826 | Detachable Rasp - Size 6 |
| 6827 | Detachable Rasp - Size 7 |
| 6828 | Detachable Rasp - Size 8 |
| 6829 | Detachable Rasp - Size 9 |
| 6830 | Detachable Rasp - Size 10 |
| 6831 | Detachable Rasp - Size 11 |
| 6832 | Detachable Rasp - Size 12 |

| | |
|---------------|--|
| 9666-22-000 | Head Trial - Size 22mm/Neutral (12/14 Taper) |
| 9666-22-350 | Head Trial - Size 22mm/+3.5mm (12/14 Taper) |
| 9666-22-800 | Head Trial - Size 22mm/+8mm (12/14 Taper) |
| 9666-26-000 | Head Trial - Size 28mm/Neutral (12/14 Taper) |
| 9666-28-004 | Head Trial - Size 28mm/-4mm (12/14 Taper) |
| 9666-28-400 | Head Trial - Size 28mm/+4mm (12/14 Taper) |
| 9666-28-800 | Head Trial - Size 28mm/+8mm (12/14 Taper) |
| 9666-32-000 | Head Trial - Size 32mm/Neutral (12/14 Taper) |
| 9666-32-004 | Head Trial - Size 32mm/-4mm (12/14 Taper) |
| 9666-32-400 | Head Trial - Size 32mm/+4mm (12/14 Taper) |
| 9666-32-800 | Head Trial - Size 32mm/+8mm (12/14 Taper) |
| 6852 | Head Trial for Detachable Rasp - 28 - Small |
| 6853 | Head Trial for Detachable Rasp - 28 - Med |
| 6854 | Head Trial for Detachable Rasp - 28 - Large |
| 01.00129-284 | Head Trial for Detachable Rasp - 28 - XL |
| 6856 | Head Trial for Detachable Rasp - 32 - Small |
| 6857 | Head Trial for Detachable Rasp - 32 - Med |
| 6858 | Head Trial for Detachable Rasp - 32 - Large |
| 01.000129.324 | Head Trial for Detachable Rasp - 32 - XL |

ALLOCLASSIC SL HIP INSTRUMENT CASE 2

| Catalog No. | Description |
|--------------|--------------------------------------|
| 01.00129.500 | Instrument Case (Alloclassic Trials) |
| 01.00029.031 | Instrument Case Cover |
| 6862 | Alloclassic SL Trial Stem, Size 01 |
| 6863 | Alloclassic SL Trial Stem, Size 0 |
| 6864 | Alloclassic SL Trial Stem, Size 1 |
| 6865 | Alloclassic SL Trial Stem, Size 2 |
| 6866 | Alloclassic SL Trial Stem, Size 3 |
| 6867 | Alloclassic SL Trial Stem, Size 4 |
| 6868 | Alloclassic SL Trial Stem, Size 5 |
| 6869 | Alloclassic SL Trial Stem, Size 6 |
| 6870 | Alloclassic SL Trial Stem, Size 7 |
| 6871 | Alloclassic SL Trial Stem, Size 8 |
| 6872 | Alloclassic SL Trial Stem, Size 9 |
| 6873 | Alloclassic SL Trial Stem, Size 10 |
| 6874 | Alloclassic SL Trial Stem, Size 11 |
| 6875 | Alloclassic SL Trial Stem, Size 12 |

ALLOCLASSIC SLL HIP IMPLANTS

| Catalog No. | Description |
|-------------|---------------------------------------|
| 2883 | Alloclassic SLL Revision Hip, Size 3 |
| 2884 | Alloclassic SLL Revision Hip, Size 4 |
| 2885 | Alloclassic SLL Revision Hip, Size 5 |
| 2886 | Alloclassic SLL Revision Hip, Size 6 |
| 2887 | Alloclassic SLL Revision Hip, Size 7 |
| 2888 | Alloclassic SLL Revision Hip, Size 8 |
| 2889 | Alloclassic SLL Revision Hip, Size 9 |
| 2890 | Alloclassic SLL Revision Hip, Size 10 |
| 2891 | Alloclassic SLL Revision Hip, Size 11 |

ALLOCLASSIC SLL HIP INSTRUMENT CASE 1

| Catalog No. | Description |
|--------------|--|
| 8095 | Striking Platform/Handle, Alloclassic |
| 01.00209.014 | Tapered Bolt (sizes 1-12), Alloclassic |
| 5112 | Extractor Hook |
| 5671 | Gliding Hammer |
| 5672 | Impactor |
| 01.00129.650 | Chisel |
| 6835 | Adapter with Lever |
| 8023 | Detachable Rasp/Revision - Size 3 |
| 8024 | Detachable Rasp/Revision - Size 4 |
| 8025 | Detachable Rasp/Revision - Size 5 |
| 8026 | Detachable Rasp/Revision - Size 6 |
| 8027 | Detachable Rasp/Revision - Size 7 |
| 8028 | Detachable Rasp/Revision - Size 8 |
| 8029 | Detachable Rasp/Revision - Size 9 |
| 8030 | Detachable Rasp/Revision - Size 10 |

| | |
|--------------|--|
| 8031 | Detachable Rasp/Revision - Size 11 |
| 9666-22-000 | Head Trial - Size 22mm/Neutral (12/14 Taper) |
| 9666-22-350 | Head Trial - Size 22mm/+3.5mm (12/14 Taper) |
| 9666-22-800 | Head Trial - Size 22mm/+8mm (12/14 Taper) |
| 9666-26-000 | Head Trial - Size 28mm/Neutral (12/14 Taper) |
| 9666-28-004 | Head Trial - Size 28mm/-4mm (12/14 Taper) |
| 9666-28-400 | Head Trial - Size 28mm/+4mm (12/14 Taper) |
| 9666-28-800 | Head Trial - Size 28mm/+8mm (12/14 Taper) |
| 9666-32-000 | Head Trial - Size 32mm/Neutral (12/14 Taper) |
| 9666-32-004 | Head Trial - Size 32mm/-4mm (12/14 Taper) |
| 9666-32-400 | Head Trial - Size 32mm/+4mm (12/14 Taper) |
| 9666-32-800 | Head Trial - Size 32mm/+8mm (12/14 Taper) |
| 6852 | Head Trial for Detachable Rasp - 28 - Small |
| 6853 | Head Trial for Detachable Rasp - 28 - Med |
| 6854 | Head Trial for Detachable Rasp - 28 - Large |
| 01.00129-284 | Head Trial for Detachable Rasp - 28 - XL |
| 6856 | Head Trial for Detachable Rasp - 32 - Small |
| 6857 | Head Trial for Detachable Rasp - 32 - Med |
| 6858 | Head Trial for Detachable Rasp - 32 - Large |
| 01.00129.324 | Head Trial for Detachable Rasp - 32 - XL |

ALLOCLASSIC TEMPLATES

| Catalog No. | Description |
|--------------|------------------------------------|
| 06.01045.000 | Templates, Alloclassic SL Hip 18% |
| 06.01084.000 | Templates, Alloclassic SLL Hip 18% |



Please refer to package inserts for complete product information, including contraindications, warnings, precautions, and adverse effects.

Contact your Zimmer Representative or visit us at www.zimmer.com.